

m|devices®

## AN080000NS - AN084001NS Nebuliser Kits

### INSTRUCTION FOR USE

#### DESCRIPTION

m|devices Nebuliser Kits are available in 4 different mask sizes; Paediatric Standard Shape, Paediatric Elongated Shape, Adolescent Standard shape and Adult Elongated Shape. An option is available with mouthpiece only.

The prescribed liquid medicine is poured into the nebuliser jar and the tubing is connected to an oxygen source. The oxygen is applied through the tubing into the nebuliser to break up medication solutions into aerosol droplets suitable for inhalation.

All kits include a 10ml jar with volume markings and are available with two oxygen tubing lengths: 2.1m and 3m.

**FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.**

#### CAUTION

- Single use only.
- DO NOT re-sterilise.
- DO NOT store at extreme temperatures and humidity, avoid direct sunlight. Handle with care.
- NON STERILE, DO NOT use if the package or product has been damaged or contaminated.
- EU Notice: any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and /or patient is established.

#### STEPS

1. Remove product from packaging and connect the tubing to the nebuliser unit.
2. Fill the nebuliser jar with the prescribed medication and carefully secure the lid over the nebuliser jar.
3. Attach the tubing and mask to the nebuliser jar.
4. Set the flow as prescribed and check for a aerosol mist.
5. Place the mask over the patient's nose and mouth with the elastic strap below the ears and around the neck.
6. Adjust the metal nose clip over the patient's nose for added comfort for the patient and to ensure a well-fitting seal.
7. The patient should continue to breathe normally until the solution has vaporised or as clinically indicated.

m|devices

27 Llewellyn Avenue  
Villawood NSW 2163  
Australia  
T +61(0)2 8718 2800  
F +61(0)2 8718 2801

 AIRWAY  
MANAGEMENT Nebuliser Kits

 **Multigate Co., Ltd.**  
7 Lanjiang Road 312000  
Shaoxing China

 **RFX+CARE International A/S**  
Bakkegaardsvej 408  
3050 Humlebaek Denmark

