

m|devices®

IR010005 Hysteroscopy Irrigation Set

INSTRUCTION FOR USE

DESCRIPTION

m|devices Hysteroscopy Irrigation Set is designed for hysteroscopy procedures in providing irrigation to the surgical cavity. Single spike configuration with a manual hand pump to enable the clinician to push fluid through to clear the surgical field to assist with visualisation during the procedure.

FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.

CAUTION

- Single use.
- DO NOT re-sterilise.
- DO NOT store at extreme temperatures and humidity, avoid direct sunlight. Handle with care.
- STERILE (EO), DO NOT use if the package or product has been damaged or contaminated.
- EU Notice: any serious incident that has occurred in relation to the device should be reported to the manufacture and the competent authority of the Member State in which the user and /or patient is established.

STEPS

1. Close click clamp.
2. Hang fluid bag and remove protective cover from bag spike outlet.
3. Remove protective cap from spike of the irrigation set and insert into fluid bag using aseptic technique.
4. To prime the irrigation line, invert flexible pump chamber and open the click clamp until the fluid flows out through the silicone tubing adaptor, checking that all air has been expelled.
5. Close click clamp.
6. Connect and secure by attaching the silicone tubing adaptor to the inline port on the endoscope/instrumentation sheath.
7. Control flow rate using the click clamp until desired flow rate is achieved.
8. Squeezing the flexible pump chamber will provide a push of fluid to flush the irrigation set if required to assist.

PRECAUTION

- Use of the flexible pump chamber should only be performed by a qualified clinician, as incorrect use could potentially cause high pressures within the surgical cavity.

m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia
T +61(0)2 8718 2800
F +61(0)2 8718 2801

 SURGICAL &
SPECIALTIES Irrigation Sets

 **Multigate Co., Ltd.**
7 Lanjiang Road 312000
Shaoxing China

 **RFX+CARE International A/S**
Bakkegaardsvej 408
3050 Humlebaek Denmark

         0123