

2024 Product Catalog

Cardinal Health™ Protexis™ Surgical Gloves

Protect what matters



 Cardinal Health™
Protexis™

Contents

3	Introduction	
4-6	Designing, manufacturing and collaborating to protect	
7	Cardinal Health Clinical Advisor Team	
8	A closer look at sustainability	
10-26	Protexis™ Surgical Glove Portfolio	
	Non-latex polyisoprene	Non-latex neoprene
10	Protexis™ PI	18 Protexis™ Neoprene
11	Protexis™ PI with Neu-Thera™	19 Protexis™ Neoprene Essential
12	Protexis™ PI Classic	
13	Protexis™ PI Blue with Neu-Thera™	Latex
14	Protexis™ PI Micro	20 Protexis™ Latex
15	Protexis™ PI Orthopaedic	21 Protexis™ Latex with Neu-Thera™
16	Protexis™ PI Ortho	22 Protexis™ Latex Classic
17	Protexis™ PI Textured	23 Protexis™ Latex Hydrogel
		24 Protexis™ Latex Blue with Neu-Thera™
		25 Protexis™ Latex Micro
		26 Protexis™ Latex Ortho
27-31	Testing standards and data	
32-37	Product quick reference	
38-39	Clinical application matrix	

Cardinal Health™ Protexis™ Surgical Gloves

From our hands to yours

As a clinician, you are entrusted with the lives and well-being of your patients every day. Cardinal Health™ Protexis™ Surgical Gloves are designed with your needs in mind, providing high-quality hand protection that you can depend on. From intentional design to manufacturing to selection, the Protexis™ Portfolio provides a variety of options and the personal support you need to help safely deliver positive clinical outcomes. Rely on Protexis™ Surgical Gloves to help protect what matters most: you and your patients.



Designed to protect



Manufactured to protect



Collaborating to protect





Designed
to protect

Designed with you in mind

Your surgical gloves need to protect you and your patients when you are performing surgery. Cardinal Health™ Protexis™ Surgical Gloves are designed with the right combination of features to support your performance and safety in the operating room. The Protexis™ Surgical Glove Portfolio offers a comprehensive selection of sizes, textures, materials and thicknesses to meet your preferences and procedural needs, so you can perform at your best.

- 1 Proprietary hand mold**
Allows for natural movement to reduce hand fatigue
- 2 Independent thumb design**
Requires minimal flexion and extension force across the palm
- 3 Neu-Thera™ Coating**
Contains glycerin, gluconolactone and provitamin B5, promoting skin moisturization*
- 4 Interlocking beaded cuff**
Helps prevent cuff roll-down during procedures





Manufactured
to protect

Protecting you, our staff and the environment

Maintaining a reliable supply of surgical gloves can be challenging, especially when you have to consider quality, sustainability and safety. That's why Cardinal Health continually invests in self-manufacturing production, capacity and sustainability — so you can simplify your supply chain and put your focus back on your patients.

85%

of clinicians preferred Protexis™ Surgical Gloves because of the quality of our product.*

50+ years

of manufacturing experience, 25 years of manufacturing excellence in Rayong, Thailand



Robust testing protocols

We monitor over 2,000 quality-control variables during the production process to help ensure each glove meets our strict quality expectations.



Investment in sustainability

Innovative manufacturing processes align with your organization's sustainability goals, helping protect our environment the same way we protect you.



Award-winning labor standards

Over a decade's worth of recognition and awards from the government of Thailand, because our staff deserve the same protection you do.



Collaborating
to protect

Our best practice is supporting you


Choosing the right surgical gloves can be difficult. The Cardinal Health™ Protexis™ Surgical Gloves Team is here to provide tailored clinical support every step of the way. We'll provide guidance in choosing the right glove for the job at hand and for your unique needs while continually equipping you with evidence-based education, resources and personalized on-site surgical glove support.

**Our team is here
to help with**




Clinical education

Stay on top of your performance goals with up-to-date, evidence-based resources.



Inventory management

Drive cost efficiencies with opportunities for SKU harmonization.



Hands-on, end-to-end clinical support

Collaborate with our team on-site for easier glove evaluation, conversion and ongoing assistance.

Cardinal Health **Clinical Advisor Team**

Providing you with products and services education

The Cardinal Health Clinical Advisor Team is comprised of skilled nurses across the country who have experience helping our customers' clinical staff convert to Protexis™ Surgical Gloves and provide ongoing support in the form of in-servicing, on-site assessments and more.

The team's focus is to provide education on products and services provided by Cardinal Health, with an emphasis on clinical best practices. Our Clinical Advisor Team has helped Protexis™ customers by:



Providing recommendations and solutions that can help drive efficiencies in the OR:

- Glove selection
- Product standardization
- Hand sensitivities
- How to safely don and doff



Supplying educational materials and insights on clinical best practices:

- On-site demonstrations
- White papers/case studies
- Continuing education courses



Supporting during product evaluation, implementation and post-implementation

Ready to get started?

Reach out to your local sales representative who can engage our Clinical Advisor Team.

A closer look at **sustainability**

Investing in what matters, you and the environment

Your trust in Protexis™ Surgical Gloves allows us to continually invest in sustainable manufacturing and supply chain practices to help ensure we can provide you with consistent, premium hand protection that's ethically produced.



17%

reduction in greenhouse gas (GHG) emissions per pair of Protexis™ Surgical Gloves manufactured¹



>90%

of manufacturing waste has been diverted from landfills, using waste reduction efforts such as recycling plastic, paper, rubber scrap and waste boiler ash²



100%

of paper wallets, dispenser boxes and shipping cases for Protexis™ Surgical Gloves packaging are recyclable³



¹ Data on file, Cardinal Health, Glove Emissions Intensity Reports FY19-23, updated April 2024

² Data on file, Cardinal Health, Waste Disposal Reports CY21-23, updated April 2024

³ Data on file, Cardinal Health, Recyclability of Primary, Secondary and Tertiary Packaging, updated April 2024



Cardinal Health™ Protexis™ Surgical Gloves Portfolio

Synthetic polyisoprene

- 10 Protexis™ PI
- 11 Protexis™ PI with Neu-Thera™
- 12 Protexis™ PI Classic
- 13 Protexis™ PI Blue with Neu-Thera™
- 14 Protexis™ PI Micro
- 15 Protexis™ PI Orthopaedic
- 16 Protexis™ PI Ortho
- 17 Protexis™ PI Textured

Synthetic neoprene

- 18 Protexis™ Neoprene
- 19 Protexis™ Neoprene Essential

Latex

- 20 Protexis™ Latex
- 21 Protexis™ Latex with Neu-Thera™
- 22 Protexis™ Latex Classic
- 23 Protexis™ Latex Hydrogel
- 24 Protexis™ Latex Blue with Neu-Thera™
- 25 Protexis™ Latex Micro
- 26 Protexis™ Latex Ortho



This product is not made with natural rubber latex

Protexis™ PI

- Designed to be comfortable and reliable for a variety of surgical procedures
- Our most popular glove in the U.S.
- Synthetic polyisoprene, not made with natural rubber latex

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D72PT55X	5.5	11.3 in./ 287 mm	9.1 mil/ 0.23 mm	6.7 mil/ 0.17 mm	6.7 mil/ 0.17 mm	Synthetic polyisoprene (PI)		Beaded/ Rolled	50	200
2D72PT60X	6									
2D72PT65X	6.5									
2D72PT70X	7	11.8 in./ 300 mm	9.1 mil/ 0.23 mm	6.7 mil/ 0.17 mm	6.7 mil/ 0.17 mm	Synthetic polyisoprene (PI)		Beaded/ Rolled	50	200
2D72PT75X	7.5									
2D72PT80X	8									
2D72PT85X	8.5									
2D72PT90X	9									

See Appendix page 31 for complete testing standards

See Appendix page 30 for chemotherapy agent permeation results



1. Thickness tested in accordance with ASTM D 3577



This product is not made with natural rubber latex

Protexis™ PI with Neu-Thera™

- Designed to be comfortable and reliable for a variety of surgical procedures
 - Same great engineering as our Protexis™ PI with our patented **Neu-Thera™** coating
 - Synthetic polyisoprene, not made with natural rubber latex
- Ⓞ Neu-Thera™ is a moisturizing coating that we place on the inside of Protexis™ PI with Neu-Thera™. It promotes overall skin well-being by moisturizing dry, flaky skin.

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73TE55	5.5	11.3 in./ 287 mm	9.1 mil/ 0.23 mm	9.8 mil/ 0.25 mm	6.7 mil/ 0.17 mm	Synthetic polyisoprene (PI) with Neu-Thera™ Emollient Coating		Beaded/ Rolled	50	200
2D73TE60	6									
2D73TE65	6.5									
2D73TE70	7									
2D73TE75	7.5									
2D73TE80	8									
2D73TE85	8.5									
2D73TE90	9	11.8 in./ 300 mm								

See Appendix page 31 for complete testing standards





This product is not made with natural rubber latex

Protexis™ PI Classic

- Designed to be comfortable and reliable for a variety of surgical procedures
- Thicker than Protexis™ PI Surgical Gloves
- Synthetic polyisoprene, not made with natural rubber latex

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D72PL55X	5.5	11.5 in./ 292 mm	11.2 mil/ 0.29 mm	8.3 mil/ 0.21 mm	7.1 mil/ 0.18 mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	50	200
2D72PL60X	6									
2D72PL65X	6.5									
2D72PL70X	7	12 in./ 305 mm								
2D72PL75X	7.5									
2D72PL80X	8									
2D72PL85X	8.5									
2D72PL90X	9									

See Appendix page 31 for complete testing standards





This product is not made with natural rubber latex

Protexis™ PI Blue with Neu-Thera™

- Distinct **blue** color aids in alerting wearers to perforations in the outer glove
- Ideal underglove when double-gloving
- Synthetic polyisoprene, not made with natural rubber latex

Ⓞ Neu-Thera™ is a moisturizing coating that we place on the inside of Protexis™ PI with Neu-Thera™. It promotes overall skin well-being by moisturizing dry, flaky skin.



Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73EB55	5.5	11.3 in./ 287 mm	7.9 mil/ 0.20 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Synthetic polyisoprene (PI) with Neu-Thera™ Emollient Coating		Beaded/ Rolled	50	200
2D73EB60	6									
2D73EB65	6.5	11.8 in./ 300 mm	7.9 mil/ 0.20 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Synthetic polyisoprene (PI) with Neu-Thera™ Emollient Coating		Beaded/ Rolled	50	200
2D73EB70	7									
2D73EB75	7.5									
2D73EB80	8									
2D73EB85	8.5									
2D73EB90	9									

See Appendix page 31 for complete testing standards

See Appendix page 30 for chemotherapy agent permeation results





This product is not made with natural rubber latex

Protexis™ PI Micro

- Thinnest glove in the Protexis™ Synthetic Polyisoprene Portfolio
- Ideal for vascular, ophthalmology, and laparoscopic and robotics procedures
- Heightened tactile response with a comfortable, smooth, anti-slip finish
- Thin for enhanced flexibility and tactile sensitivity
- Synthetic polyisoprene, not made with natural rubber latex

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73PM55	5.5	11.3 in./ 287 mm	7.9 mil/ 0.20 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	50	200
2D73PM60	6									
2D73PM65	6.5									
2D73PM70	7	11.8 in./ 300 mm								
2D73PM75	7.5									
2D73PM80	8									
2D73PM85	8.5									
2D73PM90	9									

See Appendix page 31 for complete testing standards





This product is not made with natural rubber latex

Protexis™ PI Orthopaedic

- Thickest glove in the Protexis™ Synthetic Polyisoprene Portfolio
- Smooth finish for tactile sensitivity
- Water-based hydrogel coating for easy donnability
- Durable for broaching and with the dexterity needed for pinning
- Rich brown color reduces glare from surgical lighting
- Synthetic polyisoprene, not made with natural rubber latex

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73HT60	6	11.5 in./ 292 mm	13.4 mil/ 0.35 mm	10.2 mil/ 0.26 mm	8.3 mil/ 0.21 mm	Synthetic polyisoprene (PI) with water-based hydrogel polymer coating		Beaded/ Rolled	40	160
2D73HT65	6.5									
2D73HT70	7									
2D73HT75	7.5	12.0 in./ 305 mm								
2D73HT80	8									
2D73HT85	8.5									
2D73HT90	9									

See Appendix page 31 for complete testing standards





This product is not made with natural rubber latex

Protexis™ PI Ortho

- Thinner than PI Orthopaedic¹
- Designed for use in orthopedic procedures such as trauma, labor and delivery, and reconstructive surgery
- An ideal outer glove in a double-gloving system complemented with a colored underglove

Cat. no.	Size	Length	Thickness ²			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73ET60	6	11.5 in./ 292 mm	12.2 mil/ 0.31 mm	8.3 mil/ 0.21 mm	7.5 mil/ 0.19 mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	40	160
2D73ET65	6.5									
2D73ET70	7									
2D73ET75	7.5									
2D73ET80	8									
2D73ET85	8.5									
2D72LS90	9	12.0 in./ 305 mm								

See Appendix page 31 for complete testing standards



1. As compared to Cardinal Health™ Protexis™ Latex Surgical Gloves
 2. In accordance with ASTM D 3577



This product is not made with natural rubber latex

Protexis™ PI Textured

- Textured finish provides exceptional grip under dry and damp conditions
- Designed with anatomic fit to help reduce hand fatigue and texture to enhance grip
- Innovative textured zones of the glove enhance tactile sensitivity in dry or damp conditions
- Water-based hydrogel coating for easy donning enhances second-skin comfort

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D72TG55	5.5	11.4 in./ 292 mm	10.6 mil/ 0.28 mm	8.3 mil/ 0.21 mm	7.1 mil/ 0.18 mm	Synthetic polyisoprene (PI)	Cream	Beaded/ Rolled	50	200
2D72TG60	6									
2D72TG65	6.5									
2D72TG70	7	11.9 in./ 305 mm								
2D72TG75	7.5									
2D72TG80	8									
2D72TG85	8.5									
2D72TG90	9									

See Appendix page 31 for complete testing standards





This product is not made with natural rubber latex

Protexis™ Neoprene

- Thin and soft for enhanced tactile response¹
- Synthetic neoprene, not made with natural rubber latex
- Nitrile coating for strength, protection and easy donning

Cat. no.	Size	Length	Thickness ²			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73DP55	5.5	11.1 in./ 279 mm	6.7 mil/ 0.17 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Synthetic neoprene with nitrile polymer coating		Beaded/ Rolled	50	200
2D73DP60	6									
2D73DP65	6.5									
2D73DP70	7	11.7 in./ 295 mm	6.7 mil/ 0.17 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Synthetic neoprene with nitrile polymer coating		Beaded/ Rolled	50	200
2D73DP75	7.5									
2D73DP80	8									
2D73DP85	8.5									
2D73DP90	9									

See Appendix page 31 for complete testing standards

See Appendix page 30 for chemotherapy agent permeation results



1. Compared to Duraprene SMT (previous generation)

2. Thickness tested in accordance with ASTM D 3577



This product is not made with natural rubber latex

Protexis™ Neoprene Essential

- Manufactured with zinc oxide as an alternative to traditional chemical accelerators
- Smooth finish for tactile sensitivity
- Nitrile coating for strength, protection and easy donning
- Utilizes a specific formulation of zinc oxide during the curing process as an alternative to the four classes of chemical accelerators¹
- Synthetic neoprene, not made with natural rubber latex

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73DS55	5.5	11.1 in./ 279 mm	6.7 mil/ 0.17 mm	≥ 5.5mil/ ≥ 0.14 mm	≥ 5.5mil/ ≥ 0.14 mm	Synthetic neoprene with nitrile polymer coating	 Light brown	Beaded/ Rolled	50	200
2D73DS60	6									
2D73DS65	6.5									
2D73DS70	7	11.7 in./ 295 mm								
2D73DS75	7.5									
2D73DS80	8									
2D73DS85	8.5									
2D73DS90	9									

See Appendix page 31 for complete testing standards

Not made with traditional chemical accelerators

See Appendix page 30 for chemotherapy agent permeation results





Natural rubber latex

Protexis™ Latex

- Versatile glove to be used in a wide variety of surgical environments¹
- Gloves brown tint and opacity reduces glare and provides protection that is less obtrusive and less noticeable

Cat. no.	Size	Length	Thickness ²			Material	Color	Cuff type	Protein content ³	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72NS55X	5.5	11.1 in./ 282 mm	9.8 mil/ 0.25 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer coating	 Light Brown	Beaded/ Rolled	Less than 50 µg/dm ³	50	200
2D72NS60X	6										
2D72NS65X	6.5										
2D72NS70X	7										
2D72NS75X	7.5										
2D72NS80X	8										
2D72NS85X	8.5										
2D72NS90X	9	11.6 in./ 295 mm									

See Appendix page 31 for complete testing standards



1. Compared to other latex gloves in the Protexis™ Portfolio
 2. Thickness tested in accordance with ASTM D 3577
 3. Protein content tested using ASTM D 5712



Natural rubber latex

Protexis™ Latex with Neu-Thera™

- Designed to be comfortable and reliable for a variety of surgical procedures
- Protexis™ Latex with **Neu-Thera™** Surgical Gloves deliver exceptional protection while promoting skin moisturization¹

Ⓢ Neu-Thera™ is a moisturizing coating/emollient that is placed on the inside of Protexis™ Latex with Neu-Thera™. It promotes overall skin well-being by moisturizing dry, flaky skin.

Cat. no.	Size	Length	Thickness ²			Material	Color	Cuff type	Protein content ³	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D73TP55	5.5	11.1 in./ 282 mm	9.8 mil/ 0.25 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Light brown	Beaded/ Rolled	Less than 50 µg/dm ³	50	200
2D73TP60	6										
2D73TP65	6.5										
2D73TP70	7	11.7 in./ 295 mm	9.8 mil/ 0.25 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Light brown	Beaded/ Rolled	Less than 50 µg/dm ³	50	200
2D73TP75	7.5										
2D73TP80	8										
2D73TP85	8.5										
2D73TP90	9										

See Appendix page 31 for complete testing standards



1. Data on file with Cardinal Health, California Skin Research Institute Study, Project Number 03-118 5 In accordance with ASTM D 3577.
 2. Thickness tested in accordance with ASTM D 3577
 3. Protein content tested using ASTM D 5712



Natural rubber latex

Protexis™ Latex Classic

- Designed to protect in a broad range of cases
- Ideal outer glove when double-gloving, or can be worn as a stand-alone glove
- Exceptional protection, dexterity and an advanced grip

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Protein content ²	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72N55X	5.5	11.5 in./ 292 mm	9.8 mil/ 0.25 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer coating	 Cream	Beaded/ Rolled	Less than 50 µg/dm ²	50	200
2D72N60X	6										
2D72N65X	6.5										
2D72N70X	7	12 in./ 305 mm									
2D72N75X	7.5										
2D72N80X	8										
2D72N85X	8.5										
2D72N90X	9										

See Appendix page 31 for complete testing standards



1. Thickness tested in accordance with ASTM D 3577
 2. Protein content tested using ASTM D 5712



Natural rubber latex

Protexis™ Latex Hydrogel

- Balances tactile sensitivity with protection, even when double-gloving
- Water-based hydrogel coating for easy donning with wet or dry hands and enhances second-skin comfort of latex

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Protein content ²	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72LS55	5.5	11.5 in./ 292 mm	9.8 mil/ 0.29 mm	9.1 mil/ 0.14 mm	7.1 mil/ 0.14 mm	Natural rubber latex coated with acrylic hydrogel polymer coating	Translucent Yellow	Beaded/ Rolled	Less than 50 µg/dm ²	50	200
2D72LS60	6										
2D72LS65	6.5										
2D72LS70	7	12 in./ 305 mm	9.8 mil/ 0.29 mm	9.1 mil/ 0.14 mm	7.1 mil/ 0.14 mm	Natural rubber latex coated with acrylic hydrogel polymer coating	Translucent Yellow	Beaded/ Rolled	Less than 50 µg/dm ²	50	200
2D72LS75	7.5										
2D72LS80	8										
2D72LS85	8.5										
2D72LS90	9										

See Appendix page 31 for complete testing standards





Natural rubber latex


Protexis™ Latex Blue with Neu-Thera™

- Distinct **blue** color aids in alerting wearers to perforations in the outerglove
- Ideal underglove when double-gloving
- Designed to be comfortable and reliable for a variety of surgical procedures

Ⓞ Neu-Thera™ is a moisturizing coating that we place on the inside of Protexis™ PI with Neu-Thera™. It promotes overall skin well-being by moisturizing dry, flaky skin.

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Protein content ²	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72LU55	5.5	11.1 in./ 282 mm	7.87 mil/ 0.19 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Blue	Beaded/ Rolled	Less than 50 µg/dm ²	50	200
2D72LU60	6										
2D72LU65	6.5										
2D72LU70	7	11.6 in./ 295 mm	7.87 mil/ 0.19 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Blue	Beaded/ Rolled	Less than 50 µg/dm ²	50	200
2D72LU75	7.5										
2D72LU80	8										
2D72LU85	8.5										
2D72LU90	9										

See Appendix page 31 for complete testing standards





Natural rubber latex

Protexis™ Latex Micro

- Thinnest glove in the Protexis™ Latex Portfolio¹
- Ideal for vascular, ophthalmology, and laparoscopic and robotics procedures
- Ideal in a thin double-gloving system where fingertip sensation is essential
- Protexis™ Latex Micro is approximately 20 percent thinner than Protexis™ Latex for enhanced flexibility and tactile sensitivity
- Heightened tactile response with a comfortable, smooth, anti-slip finish

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Protein content ²	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72NT55X	5.5	11.1 in./ 282 mm	6.9 mil/ 0.17 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Natural rubber latex with nitrile polymer coating	 Light brown	Beaded/ Rolled	Less than 50 µg/dm ²	50	200
2D72NT60X	6										
2D72NT65X	6.5										
2D72NT70X	7	11.6 in./ 295 mm	6.9 mil/ 0.17 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Natural rubber latex with nitrile polymer coating	 Light brown	Beaded/ Rolled	Less than 50 µg/dm ²	50	200
2D72NT75X	7.5										
2D72NT80X	8										
2D72NT85X	8.5										
2D72NT90X	9										

See Appendix page 31 for complete testing standards



1. Thickness tested in accordance with ASTM D 3577
 2. Protein content tested using ASTM D 5712



Natural rubber latex

Protexis™ Latex Ortho

- Thickest glove in the Protexis™ Latex Portfolio
- Smooth finish for tactile sensitivity
- Water-based hydrogel coating for easy donning
- Durable for broaching and tactile for pinning
- Rich brown color reduces glare from surgical lighting

Cat. no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Protein content ²	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72LT60	6	11.1 in./ 282 mm	13.4 mil/ 0.35 mm	9.4 mil/ 0.24 mm	8.3 mil/ 0.21 mm	Natural rubber latex with water-based hydrogel polymer coating	 Brown	Less than 50 µg/dm ²	40	160	
2D72LT65	6.5										
2D72LT70	7	11.6 in./ 295 mm									
2D72LT75	7.5										
2D72LT80	8										
2D72LT85	8.5										
2D72LT90	9										

See Appendix page 31 for complete testing standards



1. Thickness tested in accordance with ASTM D 3577
 2. Protein content tested using ASTM D 5712



Appendix: Cardinal Health™ Protexis™ Testing standards and technical data

Non-latex polyisoprene



This product is not made with natural rubber latex



Properties (before aging)

	Protexis™ PI	Protexis™ PI with Neu-Thera™	Protexis™ PI Classic	Protexis™ PI Blue with Neu-Thera™	Protexis™ PI Micro
Tensile strength (min)	≥ 17 MPa ¹	≥ 17 MPa ¹	≥ 17 MPa ¹	≥ 17 MPa ¹	≥ 17 MPa ¹
Stress at 500% elongation (modulus) (max)	≤ 7.0 MPa ¹	≤ 7.0 MPa ¹	≤ 7.0 MPa ¹	≤ 7.0 MPa ¹	≤ 7.0 MPa ¹
Ultimate elongation (elasticity) (min)	≥ 650% ¹	≥ 650% ¹	≥ 650% ¹	≥ 650% ¹	≥ 650% ¹
Puncture resistance (cuff) ²	AV ≥ 5N	AV ≥ 5N	AV ≥ 5N	AV ≥ 5N	AV ≥ 5N
Freedom from holes ³	0.65 AQL ¹	0.65 AQL ¹	0.65 AQL ¹	0.65 AQL ¹	0.65 AQL ¹
Sterilization	Radiation	Radiation	Radiation	Radiation	Radiation
Accelerant	Zinc diethyldithiocarbamate (ZDEC), Zinc mercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG)				



Properties (before aging)

	Protexis™ PI Orthopaedic	Protexis™ PI Ortho	Protexis™ PI Textured
Tensile strength (min)	≥ 17 Mpa ¹	≥ 17 MPa ¹	≥ 17 MPa ¹
Stress at 500% elongation (modulus) (max)	≤ 7.0 Mpa ¹	≤ 7.0 MPa ¹	≤ 7.0 MPa ¹
Ultimate elongation (elasticity) (min)	≥ 650% ¹	≥ 650% ¹	≥ 650% ¹
Puncture resistance (cuff) ²	AV = 5N	AV ≥ 5	N/A
Freedom from holes ³	0.65 AQL ¹	0.65 AQL ¹	0.65 AQL ¹
Sterilization	Radiation	Radiation	Radiation
Accelerant	Zinc diethyldithiocarbamate (ZDEC), Zinc mercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG)		

1. In accordance with ASTM D 3577


2. Tested in accordance with AS/NZS 4179, min 5 N

3. Tested in accordance with ASTM D 5151

Non-latex neoprene



This product is not made with natural rubber latex



Properties (before aging)

	Protexis™ Neoprene	Protexis™ Neoprene Essential
Tensile strength (min)	≥ 17 MPa ¹	≥ 17 MPa ¹
Stress at 500% elongation (modulus) (max)	≤ 7.0 MPa ¹	≤ 7.0 MPa ¹
Ultimate elongation (elasticity) (min)	≥ 650% ¹	≥ 650% ¹
Puncture resistance (cuff) ²	AV ≥ 5	AV ≥ 5
Freedom from holes ³	0.65 AQL ¹	0.65 AQL ¹
Sterilization	Radiation	Radiation
Accelerant	ZDBC (Zinc Dibutyldithiocarbamate)	Zinc Oxide

1. In accordance with ASTM D 3577
2. Tested in accordance with AS/NZS 4179, min 5 N
3. Tested in accordance with ASTM D 5151

Latex



NATURAL RUBBER LATEX



Properties (before aging)

	Protexis™ Latex	Protexis™ Latex with Neu-Thera™	Protexis™ Latex Classic	Protexis™ Latex Hydrogel	Protexis™ Latex Blue with Neu-Thera™	Protexis™ Latex Micro	Protexis™ Latex Ortho
Tensile strength (min)	≥ 24 MPa ¹	≥ 24 MPa ¹	≥ 24 MPa ¹	≥ 24 MPa ¹	≥ 24 MPa ¹	≥ 24 MPa ¹	≥ 24 Mpa ¹
Stress at 500% elongation (modulus) (max)	≤ 5.5 MPa ¹	≤ 7.0 MPa ¹	≤ 5.5 MPa ¹	≤ 5.5 MPa ¹	≤ 5.5 MPa ¹	≤ 5.5 MPa ¹	≤ 5.5 Mpa ¹
Ultimate elongation (elasticity) (min)	>/=750% ¹	≥ 650% ¹	≥ 750% ¹	≥ 750% ¹	≥ 750% ¹	≥ 750% ¹	≥ 750% ¹
Puncture resistance (cuff) ²	AV ≥ 5	AV ≥ 5	AV ≥ 5	AV ≥ 5	AV ≥ 5	AV ≥ 5	AV ≥ 5
Freedom from holes ³	0.65 AQL ¹	0.65 AQL ¹	0.65 AQL ¹	0.65 AQL ¹	0.65 AQL ¹	0.65 AQL ¹	0.65 AQL ¹
Sterilization	Radiation	Radiation	Radiation	Radiation	Radiation	Radiation	Radiation
Accelerant	ZDBC (Zinc Dibutyldithiocarbamate)						

1. In accordance with ASTM D 3577
2. Tested in accordance with AS/NZS 4179, min 5 N
3. Tested in accordance with ASTM D 5151

Chemotherapy agent permeation testing

All gloves listed below meet the requirements described in USP <800>, including being powder-free and meeting ASTM Standard D6978 testing for chemotherapy drug permeation.

Agent	Minimum breakthrough detection time in minutes (0.01 µg/cm ² /minute)			
	Protexis™ PI ¹	Protexis™ Neoprene ²	Protexis™ Neoprene Essential ²	Protexis™ PI Blue with Neu-Thera™ ¹
Carmustine (3.3 mg/mL)	15.26	28.2	24.6	18.5
Cisplatin (1.0 mg/mL)	> 240	> 240	> 240	> 240
Cyclophosphamide (20 mg/mL)	> 240	> 240	> 240	> 240
Doxorubicin HCL (2.0 mg/mL)	> 240	> 240	> 240	> 240
Etoposide (20 mg/mL)	> 240	> 240	> 240	> 240
Fluorouracil (50 mg/mL)	> 240	> 240	> 240	> 240
Ifosfamide (50 mg/mL)	> 240	Not tested	Not tested	Not tested
Methotrexate (25 mg/mL)	> 240	> 240	> 240	> 240
Mitomycin C (0.5 mg/ml)	> 240	Not tested	> 240	> 240
Mitoxantrone (2 mg/mL)	> 240	Not tested	Not tested	Not tested
Paclitaxel (6.0 mg/mL)	> 240	> 240	> 240	> 240
ThioTEPA (10 mg/mL)	16.04	48.9	83.1	24.4
Vincristine Sulfate (1.0 mg/mL)	> 240	> 240	> 240	> 240

Permeation times differ for gloves sterilized using gamma radiation

1. Warning: Do not use PROTEXIS™ PI or PI Blue with Ne-Thera with Carmustine (BCNU) (3.3 mg/mL) or ThioTEPA (10 mg/mL).
2. Warning: Do not use PROTEXIS™ Neoprene or Neoprene Essential with Carmustine (BCNU) (3.3 mg/mL).



When chemotherapy drugs are present, glove selection should be based on the specific type(s) of chemicals used. Users should review drug labeling or Material Safety Data Sheets for the chemicals being used to determine an adequate level of protection.

These gloves have been tested for resistance to permeation of various chemotherapy drugs per ASTM D 6978, "Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs."

Testing standards

Global Quality Standards tested for and adhered to (results on file)

ASTM D3577, EN 455-2, ISO 10282	Physical dimension (length, width, palm)
ASTM D3577, EN 455-2, ISO 10282	Physical properties (tensile strength)
ASTM D624	Tear strength testing (T-tear, V-tear)
AS/NZA 4179	Puncture resistant testing
ASTM D5151, EN 455-1, ISO 10282	Freedom from holes (water-tightness)
ASTM D6124, EN ISO 21171	Powder residue for powder-free gloves
ASTM D6124	Powder amount for powdered gloves
ASTM D5712	Aqueous extractable protein content
ASTM D6499, ASTM D3577	Antigenic protein content
EN 455-3	Leachable protein level, modified Lowry method
ASTM D7102, EN 455-3	Endotoxin
ASTM D7160	Storage stability, accelerated aging
EN 374-5/ISO 16604	Storage stability, real-time aging
ASTM F739, EN 16523	Lab chemical permeation
ASTM D6978, EN 16523	Chemotherapy drug permeation (results on page 28)
EN 374-5, ISO 16604	Bacteriophage penetration
ISO 10993-10	Sensitivity testing and primary skin irritation
ISO 10993-5	Cytotoxicity testing
ISO 10993-11	Systemic Toxicity Testing (Acute)
ASTM D3577	Sterility test

The standards only apply as relevant to the respective product family.

PPE required testing (results on file)

EN ISO 21420	Protective glove - general requirements and test method
EN ISO 374-1:2016	Performance requirements for chemicals risk
EN 374-2:2019	Resistance to penetration against dangerous chemicals and micro-organisms
EN 16523-1:2015	Materials resistance to permeation by chemicals
EN 374-4:2019	Resistance to degradation by chemicals
EN ISO 374-5:2016	Performance requirements for micro-organisms risk

Internal procedures (results on file)

- Bone cement permeation
- Low-hydration conductivity
- Residue accelerator test



Non-latex polyisoprene



This product is not made with natural rubber latex

	Cat. no.	Size	Length	Thickness*			Material <i>Coating in red</i>	Color	Cuff type	Chemical accelerant	Qty/ bx	Qty/ cs
				Finger	Palm	Cuff						
 Protexis™ PI	2D72PT55X	5.5	11.3 in./ 287 mm	9.1 mil/ 0.23 mm	6.7 mil/ 0.17 mm	6.7 mil/ 0.17 mm	Synthetic polyisoprene (PI)		Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
	2D72PT60X	6										
	2D72PT65X	6.5										
	2D72PT70X	7	11.8 in./ 300 mm									
	2D72PT75X	7.5										
	2D72PT80X	8										
	2D72PT85X	8.5										
	2D72PT90X	9										
 Protexis™ PI with Neu-Thera™	2D73TE55	5.5	11.3 in./ 287 mm	9.1 mil/ 0.23 mm	9.8 mil/ 0.25 mm	6.7 mil/ 0.17 mm	Synthetic polyisoprene (PI) with Neu-Thera™ Emollient Coating		Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
	2D73TE60	6										
	2D73TE65	6.5										
	2D73TE70	7	11.8 in./ 300 mm									
	2D73TE75	7.5										
	2D73TE80	8										
	2D73TE85	8.5										
	2D73TE90	9										
 Protexis™ PI Classic	2D72PL55X	5.5	11.5 in./ 292 mm	11.2 mil/ 0.29 mm	8.3 mil/ 0.21 mm	7.1 mil/ 0.18 mm	Synthetic polyisoprene (PI)		Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
	2D72PL60X	6										
	2D72PL65X	6.5										
	2D72PL70X	7	12 in./ 305 mm									
	2D72PL75X	7.5										
	2D72PL80X	8										
	2D72PL85X	8.5										
	2D72PL90X	9										

*Thickness tested in accordance with ASTM D 3577



This product is not made with natural rubber latex

Non-latex polyisoprene

	Cat. no.	Size	Length	Thickness*			Material <i>Coating in red</i>	Color	Cuff type	Chemical accelerant	Qty/ bx	Qty/ cs
				Finger	Palm	Cuff						
 Protexis™ PI Blue with Neu-Thera™	2D73EB55	5.5	11.3 in./ 287 mm	7.9 mil/ 0.20 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Synthetic polyisoprene (PI) with Neu-Thera™ Emollient Coating	 Blue	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
	2D73EB60	6										
	2D73EB65	6.5										
	2D73EB70	7	11.8 in./ 300 mm									
	2D73EB75	7.5										
	2D73EB80	8										
	2D73EB85	8.5										
2D73EB90	9											
 Protexis™ PI Micro	2D73PM55	5.5	11.3 in./ 287 mm	7.9 mil/ 0.20 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
	2D73PM60	6										
	2D73PM65	6.5										
	2D73PM70	7	11.8 in./ 300 mm									
	2D73PM75	7.5										
	2D73PM80	8										
	2D73PM85	8.5										
2D73PM90	9											
 Protexis™ PI Orthopaedic	2D73HT60	6	11.5 in./ 292 mm	13.4 mil/ 0.35 mm	10.2 mil/ 0.26 mm	8.3 mil/ 0.21 mm	Synthetic polyisoprene (PI) with water-based hydrogel polymer coating	 Brown	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	40	160
	2D73HT65	6.5										
	2D73HT70	7	12.0 in./ 305 mm									
	2D73HT75	7.5										
	2D73HT80	8										
	2D73HT85	8.5										
2D73HT90	9											

*Thickness tested in accordance with ASTM D 3577

Non-latex polyisoprene



This product is not made with natural rubber latex

	Cat. no.	Size	Length	Thickness*			Material <i>Coating in red</i>	Color	Cuff type	Chemical accelerant	Qty/ bx	Qty/ cs
				Finger	Palm	Cuff						
 Protexis™ PI Ortho	2D73ET60	6	11.5 in./ 292 mm	12.2 mil/ 0.31 mm	8.3 mil/ 0.21 mm	7.5 mil/ 0.19 mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	40	160
	2D73ET65	6.5										
	2D73ET70	7	12.0 in./ 305 mm									
	2D73ET75	7.5										
	2D73ET80	8										
	2D73ET85	8.5										
	2D72LS90	9										
 Protexis™ PI Textured	2D72TG55	5.5	11.4 in./ 292 mm	10.6 mil/ 0.28 mm	8.3 mil/ 0.21 mm	7.1 mil/ 0.18 mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
	2D72TG60	6										
	2D72TG65	6.5	11.9 in./ 305 mm									
	2D72TG70	7										
	2D72TG75	7.5										
	2D72TG80	8										
	2D72TG85	8.5										
	2D72TG90	9										

*Thickness tested in accordance with ASTM D 3577

Non-latex neoprene



This product is not made with natural rubber latex

	Cat. no.	Size	Length	Thickness*			Material <i>Coating in red</i>	Color	Cuff type	Chemical accelerant	Qty/ bx	Qty/ cs
				Finger	Palm	Cuff						
 Protexis™ Neoprene	2D73DP55	5.5	11.1 in./ 279 mm	6.7 mil/ 0.17 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Synthetic neoprene with nitrile polymer	 Light Brown	Beaded/ Rolled	Zinc Dibutyldithiocarbamate (ZDBC)	50	200
	2D73DP60	6										
	2D73DP65	6.5										
	2D73DP70	7	11.7 in./ 295 mm									
	2D73DP75	7.5										
	2D73DP80	8										
	2D73DP85	8.5										
2D73DP90	9											
 Protexis™ Neoprene Essential	2D73DS55	5.5	11.1 in./ 279 mm	6.7 mil/ 0.17 mm	≥ 5.5mil/ ≥ 0.14 mm	≥ 5.5mil/ ≥ 0.14 mm	Synthetic neoprene with nitrile polymer	 Light Brown	Beaded/ Rolled	Manufactured with Zinc Oxide, an alternative to traditional chemical accelerators	50	200
	2D73DS60	6										
	2D73DS65	6.5										
	2D73DS70	7										
	2D73DS75	7.5										
	2D73DS80	8										
	2D73DS85	8.5										
	2D73DS90	9										

*Thickness tested in accordance with ASTM D 3577

Latex



	Cat. no.	Size	Length	Thickness*			Material <i>Coating in red</i>	Color	Cuff type	Protein content	Chemical accelerant	Qty/ bx	Qty/ cs
				Finger	Palm	Cuff							
 Protexis™ Latex	2D72NS55X	5.5	11.1 in./ 282 mm	9.8 mil/ 0.25 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer	 Light Brown	Beaded/ Rolled	Less than 50 µg/dm ²	Zinc Dibutyldithio- carbamate (ZDBC)	50	200
	2D72NS60X	6											
	2D72NS65X	6.5											
	2D72NS70X	7	11.6 in./ 295 mm										
	2D72NS75X	7.5											
	2D72NS80X	8											
	2D72NS85X	8.5											
	2D72NS90X	9											
 Protexis™ Latex with Neu-Thera™	2D73TP55	5.5	11.1 in./ 282 mm	9.3 mil/ 0.24 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Light Brown	Beaded/ Rolled	Less than 50 µg/dm ²	Zinc Dibutyldithio- carbamate (ZDBC)	50	200
	2D73TP60	6											
	2D73TP65	6.5											
	2D73TP70	7	11.7 in./ 295 mm										
	2D73TP75	7.5											
	2D73TP80	8											
	2D73TP85	8.5											
	2D73TP90	9											
 Protexis™ Latex Classic	2D72N55X	5.5	11.5 in./ 292 mm	9.8 mil/ 0.25 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer	 Cream	Beaded/ Rolled	Less than 50 µg/dm ²	Zinc Dibutyldithio- carbamate (ZDBC)	50	200
	2D72N60X	6											
	2D72N65X	6.5											
	2D72N70X	7	12 in./ 305 mm										
	2D72N75X	7.5											
	2D72N80X	8											
	2D72N85X	8.5											
	2D72N90X	9											
 Protexis™ Latex Hydrogel	2D72LS55	5.5	11.5 in./ 292 mm	9.8 mil/ 0.29 mm	9.1 mil/ 0.14 mm	7.1 mil/ 0.14 mm	Natural rubber latex with acrylic hydrogel polymer	 Translucent Yellow	Beaded/ Rolled	Less than 50 µg/dm ²	Zinc Dibutyldithio- carbamate (ZDBC)	50	200
	2D72LS60	6											
	2D72LS65	6.5											
	2D72LS70	7	12 in./ 305 mm										
	2D72LS75	7.5											
	2D72LS80	8											
	2D72LS85	8.5											
	2D72LS90	9											

*Thickness tested in accordance with ASTM D 3577

Latex



	Cat. no.	Size	Length	Thickness*			Material <i>Coating in red</i>	Color	Cuff type	Protein content	Chemical accelerant	Qty/ bx	Qty/ cs
				Finger	Palm	Cuff							
 Protexis™ Latex Blue with Neu-Thera™	2D72LU55	5.5	11.1 in./ 282 mm	7.87 mil/ 0.19 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Blue	Beaded/ Rolled	Less than 50 µg/dm ²	Zinc Dibutylthio-carbamate (ZDBC)	50	200
	2D72LU60	6											
	2D72LU65	6.5											
	2D72LU70	7	11.6 in./ 295 mm										
	2D72LU75	7.5											
	2D72LU80	8											
	2D72LU85	8.5											
	2D72LU90	9											
 Protexis™ Latex Micro	2D72NT55X	5.5	11.1 in./ 282 mm	6.9 mil/ 0.17 mm	5.5 mil/ 0.14 mm	5.5 mil/ 0.14 mm	Natural rubber latex with nitrile polymer	 Light Brown	Beaded/ Rolled	Less than 50 µg/dm ²	Zinc Dibutylthio-carbamate (ZDBC)	50	200
	2D72NT60X	6											
	2D72NT65X	6.5											
	2D72NT70X	7	11.6 in./ 295 mm										
	2D72NT75X	7.5											
	2D72NT80X	8											
	2D72NT85X	8.5											
	2D72NT90X	9											
 Protexis™ Latex Ortho	2D72LT60	6	11.1 in./ 282 mm	13.4 mil/ 0.35 mm	9.4 mil/ 0.24 mm	8.3 mil/ 0.21 mm	Natural rubber latex with water-based hydrogel polymer coating	 Brown	Beaded/ Rolled	Less than 50 µg/dm ²	Zinc Dibutylthio-carbamate (ZDBC)	40	160
	2D72LT65	6.5											
	2D72LT70	7	11.6 in./ 295 mm										
	2D72LT75	7.5											
	2D72LT80	8											
	2D72LT85	8.5											
	2D72LT90	9											

*Thickness tested in accordance with ASTM D 3577

Synthetic portfolio (Polyisoprene and neoprene)



	Department															
	General	Cardiovascular	Dental/Maxillofacial	Ear, Nose and Throat (ENT)	Endovascular	Labor & Delivery	Laparoscopic/Robotics	Neuro	Obstetrics	Ophthalmology	Orthopedics	Pediatrics	Plastics	Thoracic	Urology	Vascular
Protexis™ PI Blue with Neu-Thera™ 2D73EB55-90	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀
Protexis™ PI Micro 2D73PM55-90	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀
Protexis™ PI 2D72PT55X-90X	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀
Protexis™ PI with Neu-Thera™ 2D73TE55-90	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀
Protexis™ PI Classic 2D72PL55X-90X	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀
Protexis™ PI Orthopaedic 2D73HT60-90	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀
Protexis™ PI Ortho 2D73ET60-90	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀
Protexis™ PI Textured 2D72TG55-90	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀
Protexis™ Neoprene 2D73DP55-90	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀
Protexis™ Neoprene Essential* 2D73DS55-90	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀	▶▶◀◀

Key

Glove thickness level

thin (0.14-0.20 mm) ▶▶◀◀

standard (0.20-0.28 mm) ▶▶◀◀

thick (0.24-0.34 mm) ▶▶◀◀

Ranges from finger and palm thickness

Smooth |

Moderate ||

High |||

Highest ||||

Grip levels

Non-textured ○

Textured ●

Moisture coating

Neu-Thera™ emollient coating for hand health

Triple dip nitrile coating for strength

Hydrogel coating

For easier donning, wet or dry

Standard

Underglove

Glove style

Alternative chemical accelerator*

Chemotherapy agent permeation testing†

Chemical agent permeation testing‡

This table was developed by a group of clinicians. It reflects current best practices of surgical glove usage per application. Ultimately, it is up to the discretion of the clinician to choose the right glove for the procedure

* Uses a specific formulation of zinc oxide during the curing process as an alternative to traditional chemical accelerators

† Meets the ASTM D 6978 standard – When chemotherapy drugs are present, gloves selection should be based on the specific type(s) of chemicals used. Users are recommended to review drug labeling or material safety data sheets for the chemicals being used to determine an adequate level of protection.

‡ Meets the ASTM F739 standard – When chemical drugs are present, gloves selection should be based on the specific type(s) of chemicals used. Users are recommended to review drug labeling or material safety data sheets for the chemicals being used to determine an adequate level of protection.

Latex portfolio



Caution: This product contains natural rubber latex which may cause allergic reactions.
Caution: Safe use of these gloves by latex sensitized individuals has not been established

	Department	Department															
		General	Cardiovascular	Dental/Maxillofacial	Ear, Nose and Throat (ENT)	Endovascular	Labor & Delivery	Laparoscopic/Robotics	Neuro	Obstetrics	Ophthalmology	Orthopedics	Pediatrics	Plastics	Thoracic	Urology	Vascular
Protexis™ Latex Blue with Neu-Thera™ 2D72LU55-90	Thin, Water, Smooth, Underglove, No Chem, No Chem		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Protexis™ Latex Micro 2D72NT55X-90X	Thin, Water, Smooth, Standard, No Chem, No Chem				●	●		●				●	●				●
Protexis™ Latex 2D72NS55X-90X	Standard, Water, Smooth, Underglove, No Chem, No Chem	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Protexis™ Latex with Neu-Thera™ 2D73TP55-90	Standard, Water, Smooth, Underglove, No Chem, No Chem	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Protexis™ Latex Classic 2D72N55X-90X	Standard, Water, Smooth, Underglove, No Chem, No Chem	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
Protexis™ Latex Hydrogel 2D72LS55-90	Standard, Hydrogel, Smooth, Underglove, No Chem, No Chem	●	●	●				●	●	●	●			●	●		
Protexis™ Latex Ortho 2D72LT60-90	Thick, Water, Smooth, Underglove, No Chem, No Chem	●	●	●				●	●	●	●			●	●		

This table was developed by a group of clinicians. It reflects current best practices of surgical glove usage per application. Ultimately, it is up to the discretion of the clinician to choose the right glove for the procedure

Key

Glove thickness level

- thin (0.14-0.20 mm)
- standard (0.20-0.28 mm)
- thick (0.24-0.34 mm)

Grip levels

- Smooth
- Moderate
- High
- Highest

Texture

- Non-textured
- Textured

Moisture coating

- Neu-Thera™ emollient coating for hand health
- Triple dip nitrile coating for strength
- Hydrogel coating

For easier donning, wet or dry

- Standard
- Underglove

Alternative chemical accelerator*

Chemotherapy agent permeation testing†

Chemical agent permeation testing‡

* Uses a specific formulation of zinc oxide during the curing process as an alternative to traditional chemical accelerators

† Meets the ASTM D 6978 standard – When chemotherapy drugs are present, gloves selection should be based on the specific type(s) of chemicals used. Users are recommended to review drug labeling or material safety data sheets for the chemicals being used to determine an adequate level of protection.

‡ Meets the ASTM F739 standard – When chemical drugs are present, gloves selection should be based on the specific type(s) of chemicals used. Users are recommended to review drug labeling or material safety data sheets for the chemicals being used to determine an adequate level of protection.



Protect
what
matters



CardinalHealth™

cardinalhealth.com/protexis

© 2024 Cardinal Health. All rights reserved. CARDINAL HEALTH, the Cardinal Health LOGO, NEU-THERA and PROTEXIS are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Lit. No. 2MS24-2980625 (09/2024)