


Laryngeal Mask Airways

Our smooth one-piece design allows for easy insertion and the guide markings on the handle assist with positioning.

The translucent tubing provides clear visibility. The soft-seal cuff reduces pressure and minimises potential trauma. Our masks also feature a 15mm connector and pilot balloon to control the cuff.

- ⊗ One-piece design for easy insertion
- ⊗ Clear guide markers for easy positioning
- ⊗ Soft seal cuff to reduce pressure and minimise trauma
- ⊗ Pilot balloon to aid visual and tactile assessment of cuff
- ⊗ 100% medical grade silicone
- ⊗ Translucent tubing for clear visual monitoring
- ⊗ Available in adult, paediatric and neonate sizes
- ⊗ Colour coded for simple size selection



Laryngeal Mask Airways				STERILE
RE-ORDER	DESCRIPTION	SIZE	COLOUR	CARTON QTY
AN030000	Silicone Disposable Laryngeal Airway Mask	1	Purple	10
AN030001	Silicone Disposable Laryngeal Airway Mask	1.5	White	10
AN030002	Silicone Disposable Laryngeal Airway Mask	2	Red	10
AN030003	Silicone Disposable Laryngeal Airway Mask	2.5	Orange	10
AN030004	Silicone Disposable Laryngeal Airway Mask	3	Green	10
AN030005	Silicone Disposable Laryngeal Airway Mask	4	Yellow	10
AN030006	Silicone Disposable Laryngeal Airway Mask	5	Pink	10