

Nebuliser Kit

PRODUCT DATA SHEET



DESCRIPTION

m|devices Nebuliser Kits are available in 4 different mask sizes; Paediatric Standard Shape, Paediatric Elongated Shape, Adolescent Standard shape and Adult Elongated Shape. A mouth piece only kit is also available. They all include 2.1m and 3m kink resistant star lumen tubing and a 10mL jar. All kits are made from medical grade PVC.

FEATURES AND BENEFITS

- 4 sizes available ranging from Paediatric to Adult
- Standard shape and Elongated shape options
- 10mL jar included
- Kink resistant star lumen tubing
- Mouth piece versions available AN082000NS
- Standard connectors

PRODUCT DATA

- Tubing length: 2.1m, 3m
- Tubing type: Star Lumen
- Tubing diameter: 3.94mm
- Material: PVC
- DEHP: Yes
- Latex: No
- Shelf Life: 5 Years
- ARTG Number: 228231

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	TUBING LENGTH	SIZE	CARTON/QTY
AN080000NS	Nebuliser Kit - Standard Shape 10mL NebuliserJar	2.1m	Paediatric	50
AN080001NS	Nebuliser Kit - Elongated Shape 10mL NebuliserJar	2.1m	Paediatric	50
AN081000NS	Nebuliser Kit - Standard Shape 10mL NebuliserJar	2.1m	Adolescent	50
AN081001NS	Nebuliser Kit - Elongated Shape 10mL NebuliserJar	2.1m	Adult	50
AN082000NS	Nebuliser Kit - Mouth Piece 10mL NebuliserJar	2.1m	-	50
AN083000NS	Nebuliser Kit - Standard Shape 10mL NebuliserJar	3m	Paediatric	50
AN083001NS	Nebuliser Kit - Elongated Shape 10mL NebuliserJar	3m	Paediatric	50
AN084000NS	Nebuliser Kit - Standard Shape 10mL NebuliserJar	3m	Adolescent	50
AN084001NS	Nebuliser Kit - Elongated Shape 10mL NebuliserJar	3m	Adult	50
AN085000NS	Nebuliser Kit - Mouth Piece 10mL NebuliserJar	3m	-	50