


**Irrigation Set - Hysteroscopy**

Our Hysteroscopy Irrigation Set is simple and easy to use. Each is provided with a single spike configuration, a manual hand pump made from flexible soft material with a stepped connector and a silicone cover end to connect to the hysteroscope.

- ⊗ Length 265cm
- ⊗ Single spike configuration
- ⊗ Manual pump
- ⊗ Click clamp



Irrigation Set - Hysteroscopy			STERILE
RE-ORDER	DESCRIPTION	TOTAL LENGTH	CARTON QTY
IR010005	Hysteroscopy Irrigation Set	260cm	30