

Microbore Extension Set 75cm

DESCRIPTION

The Microbore Extension Set is a microbore tubing with female to male luer lock ends and rotating collar, which is 75cm in length.

FEATURES & BENEFITS

- Female to male luer lock
- Rotating collar
- Microbore tubing
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter 1mm
- Tubing length 75cm
- Priming volume: 0.59mL
- Sterile EO
- 290psi rated
- The product is MRI and CT suitable
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected [Result less than the Method Detection Limit EN 14372]
- Material: PVC
- Shelf Life: 3 Years
- ARTG Number: 228230

ISO STANDARDS

- EN ISO 8536-4:2013 Infusion equipment for medical use- Part 4:Infusion sets
- EN ISO 594-1 Conical fittings with a 6% [Luer] taper for syringes, needles and certain other medical equipment – Part 1: General requirements
- EN ISO 594-2 Conical fittings with 6 % [Luer] taper for syringes, needles and certain other medical equipment – Part 2:Lock fittings

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Male Luer Lock
4	Cap
5	Tubing

Microbore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV050001	Microbore Extension Set	75cm	100