

# Cystoscopy Irrigation Sets

## DESCRIPTION

The m|devices Cystoscopy Irrigation Set is designed for Cystoscopy or Bladder Irrigation. It is made from medical grade non-toxic PVC and utilises a soft flexible drip chamber and delivers high flow rates. This comes with a male luer lock (MLL) connection end to meet all clinical requirements.

## FEATURES & BENEFITS

- Soft pillow drip chamber
- Single spike
- Rigid roller regulator for accurate flow
- Male luer lock connector


## PRODUCT DATA

- Tubing inner diameter: 4.8mm
- Tubing outer diameter: 6.8mm
- Priming volume: 64mL/71mL
- Material: PVC
- Not Made with Natural Rubber Latex
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-(2-ethylhexyl)-phthalate (DEHP).
- ARTG Number: 314960
- Shelf life: 5 Years

## STANDARDS

- EN ISO 10993-5:2009 Biological evaluation of medical devices - Tests for in vitro cytotoxicity.
- EN ISO 10993-7:2008 Biological evaluation of medical devices - Ethylene oxide sterilisation residuals
- EN ISO 14971:2012 Application of Risk Management.

## PRODUCT DIAGRAM



ITEM	DESCRIPTION		
1	Cover	5	Roller body
2	Spike	6	Roller wheel
3	Tube	7	Tube
4	Soft Chamber	8	Rotating Luer Lock

## Cystoscopy Irrigation Set

RE-ORDER	DESCRIPTION	TOTAL LENGTH	CARTON/QTY
UR020002	Cystoscopy Irrigation Set with MLL	285cm	50