

FLO CARE® FEEDING SETS & TWO PACK CONNECTOR

All Flo care® Feeding Sets are ENFit™ compatible and designed to safely and securely connect enteral tube feeds to feeding tubes and accessories.

FEATURES

- Flo care® pack feeding sets connect to enteral feed via **ENPlus cross-spike end**: for increased patient safety.
- Flo care® mobile feeding sets are shorter (1580mm) than standard feeding sets (2300mm): for patient convenience and to support mobile feeding.
- The Flo care® two pack connector is a short Y-shaped connector with two **ENPlus cross-spike ends**: enables connection to two packs of enteral tube feed and/or water at the same time for patient convenience. Quick release clamp may be used to regulate flow.
- All Flo care® feeding sets are compatible with Flo care® containers (40441, 40455, 70139).
- Medication port compatible with ENFit™ syringes: for increased patient safety.
- Made from Latex-free and DEHP-free materials.

STANDARD FEEDING SETS FOR FLO CARE® INFINITY™ PUMPS

For use with ENPlus Enteral Packs & OpTri Bottles

Flo care® Infinity™ Pack Set – Y Port



Product Code: **95347**

Units Per Carton: 30

Flo care® Infinity™ Pack Set – No Port



Product Code: **95349**

Units Per Carton: 30

Flo care® Infinity™ Pack Set – 3-way Port & Drip Chamber



Product Code: **86514**

Units Per Carton: 30

For use with Enteral Bottles

Flo care® Infinity™ Bottle Set – Y Port & Drip Chamber



Product Code: **101587**

Units Per Carton: 30

Flo care® Infinity™ Bottle Set – No Port & Drip Chamber



Product Code: **90086**

Units Per Carton: 30

Flo care® Infinity™ Bottle Set – 3-way Port & Drip Chamber



Product Code: **90084**

Units Per Carton: 30

MOBILE FEEDING SETS FOR FLO CARE® INFINITY™ PUMPS

For use with ENPlus Enteral Packs & OpTri Bottles

Flocare® Infinity™ Pack Mobile Set – Y Port



Product Code: **101586**

Units Per Carton: 30

Flocare® Infinity™ Pack Mobile Set – No Port



Product Code: **86484**

Units Per Carton: 30

GRAVITY FEEDING SETS

Not for use with Enteral Feeding Pumps

Flocare® Gravity Pack Set – Y Port & Drip Chamber



Product Code: **86460**

Units Per Carton: 30


Flocare® Gravity Pack & Bottle Set – Y Port & Drip Chamber



Product Code: **89755**

Units Per Carton: 30

FLO CARE® TWO PACK CONNECTOR



Product Code: **69915**

Units Per Carton: 30

Directions for Use

- Flocare® feeding sets are intended for 'single use only' and should be changed every 24 hours.
- Flocare® two pack connectors are intended for 'single use only' and should be changed every 24 hours.

Storage

- Store in a clean, dry environment.
- Protect from dust, rain and direct sunlight.

Ordering Information

To order contact Nutricia Customer Care **1800 889 480**.

For more information call the **Nutricia Careline 1800 438 500**

NUTRICIA
LIFE-TRANSFORMING NUTRITION

Nutricia Australia Pty Ltd.
Lvl 4, Building D, 12-24 Talavera Road,
Macquarie Park, NSW, 2113.
www.nutriciamedical.com.au

NC4845-Oct20 / NUT1180C