

NeutralSite™ Valve

Our NeutralSite™ Valve helps to prevent catheter occlusions by reducing blood reflux upon disconnection, in addition to preventing catheter related blood stream infections (CRBSI). Its flat swabable pre-slit septum allows for easy disinfection. In addition, the internal straight fluid pathway helps to minimise the priming volume and address dead space concerns. The NeutralSite™ is recommended for use up to 7 days.

- ⊗ Flat swabable pre-slit septum
- ⊗ Straight internal fluid pathway
- ⊗ Low priming volume
- ⊗ Gapless space between septum and casing providing a microbial barrier against contamination
- ⊗ Lipid resistant and free of Latex and DEHP
- ⊗ Transparent housing for clear visibility
- ⊗ Compatible with luer lock and luer slip connectors



NeutralSite™ Valve

STERILE

RE-ORDER	DESCRIPTION	CARTON QTY
IV200005	NeutralSite™ Valve	100

m|devices®

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321