

# Leg Bag

## T-Tap with Elastic Straps

### DESCRIPTION

m|devices Leg Bags are designed to be unobtrusive offering discreet storage under clothing. The various volume capacity sizes available are 350mL, 500mL and 750mL bags, non-return valve and with three optional tube lengths of 10cm, 30cm and 50cm with bonded stepped connector, cater for all user requirements. The T-Tap mechanism ensures emptying the bag can be performed quickly and easily whilst the elastic straps with button hole type securement, allows an adjustable custom fit every time. The non-woven backing of the leg bag provides comfort to the wearer.

### FEATURES & BENEFITS

- Discreet underneath clothing
- Easy to use T-Tap to empty bag
- Available in 350mL, 500mL & 750mL bag sizes
- Bonded stepped connector to tubing
- Non-return valve
- Tubing lengths: 10cm, 30cm & 50cm
- Comfortable elastic straps with button hole type securement
- Soft non-woven backing


### PRODUCT DATA

- Not Made with Natural Rubber Latex
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-[2-ethylhexyl]-phthalate (DEHP).
- Material: PVC, PP, PE & Silicone
- Straps: Elastic strap with button hole - Polyester and Elastic
- Packaging: Paper and PE
- Sterile: EO
- Shelf life: 5 years
- ARTG Number: 390594

### ISO STANDARDS

- ISO 8669-2 Urine collection bags – Part 2: Determination of dimensions
- ISO 10993 Biological evaluation of medical device
- ISO 11135 Sterilisation of health care products - Ethylene oxide

### PRODUCT DIAGRAM



ITEM	DESCRIPTION	5	Outlet Tubing
1	Protective Cap	6	T-valve
2	Connector	7	Non Return Film
3	Inlet Tubing	8	Silicone Tube
4	Bag Body	9	Elastic Strap

# Leg Bag

## T-Tap with Elastic Straps

Leg Bags				STERILE
RE-ORDER	DESCRIPTION	LENGTH (CM)	COLOUR	CARTON/QTY
<b>UR235010</b>	LEG BAG - T-Tap NRV Short Tube with Elastic Straps	10cm	350mL	100
<b>UR235030</b>	LEG BAG - T-Tap NRV Long Tube with Elastic Straps	30cm	350mL	100
<b>UR235050</b>	LEG BAG - T-Tap NRV Long Tube with Elastic Straps	50cm	350mL	100
<b>UR250010</b>	LEG BAG - T-Tap NRV Short Tube with Elastic Straps	10cm	500mL	100
<b>UR250030</b>	LEG BAG - T-Tap NRV Long Tube with Elastic Straps	30cm	500mL	100
<b>UR250050</b>	LEG BAG - T-Tap NRV Long Tube with Elastic Straps	50cm	500mL	100
<b>UR275010</b>	LEG BAG - T-Tap NRV Short Tube with Elastic Straps	10cm	750mL	100
<b>UR275030</b>	LEG BAG - T-Tap NRV Long Tube with Elastic Straps	30cm	750mL	100
<b>UR275050</b>	LEG BAG - T-Tap NRV Long Tube with Elastic Straps	50cm	750mL	100