



STRONG & ABSORBENT WATERPROOF PROTECTION HALYARD* Protectapad

TRUSTED PROTECTION

FOR SURFACES & PATIENTS

SO YOU CAN CARE WITH CONFIDENCE

Looking for a product that helps to protect your patients' clothing as well as surfaces from all types of liquid such as blood, saliva, water, chemical seepage and spills? A product that's versatile, waterproof and absorbent?

OUR SOLUTION IS HALYARD* PROTECTAPAD

HALYARD* Protectapad is constructed from four highly absorbent layers of tissue and a layer of strong polyethylene backing, laminated together to help prevent strikethrough. The embossed pattern increases liquid absorbency, helping to prevent run off. Softness ensures optimum patient comfort and also assists with cushioning and protecting instruments. HALYARD* Protectapad is laminated for extra strength wet or dry and is single-use, making it the hygienic choice to protect patients and surfaces from staining and other spills.

APPLICATIONS

- Dental bib
- Hygienic surface for dental tools
- Medical procedures underlay
- Covering benches and trays
- Under IV arm of patients
- Pillow protector
- Surface liner for hairdresser/beautician products
- Patient bib
- Surface protector

DEPARTMENTS

- Dentists
- Oncology day units
- Physiotherapists
- GPs
- Optometrists
- Laboratories
- Beauty, hair & nail salons
- Pathology labs
- Aged care facilities
- Hospitals
- Blood banks



HALYARD* Protectapad

- ✓ **HIGHLY ABSORBENT TISSUE**
Absorbs fluid quickly
- ✓ **FOUR LAYERS OF SOFT TISSUE**
Provides thickness where you need it, for optimal patient comfort and to cushion and protect instruments
- ✓ **EMBOSSED PATTERN**
Increases liquid absorbency, helping to prevent run off
- ✓ **STRONG POLY BASE**
Provides good wet and dry strength¹ and helps to prevent strikethrough to keep patients and surfaces dry
- ✓ **DISPOSABLE**
Hygienic, helps reduce risk of cross-contamination

DISPENSING

HALYARD* Protectapad has a cut sheet presentation for easy, economical dispensing. And it's available in 2 sizes, small and large – so it helps better meet your needs. The small size contains 8 packs of 100 sheets, and the large size contains 4 packs of 100 sheets, with each individual pack wrapped in plastic for hygienic storage. This material is also available in a continuous roll: HALYARD* Benchroll (Code 2911).



PRODUCT INFORMATION

| Code | Description | Size | Packaging |
|------|-----------------------------|-----------------|--|
| 2706 | HALYARD* Protectapad, Small | 28.5cm x 21.5cm | 100 sheets/pack 8 packs/case 800 sheets/case |
| 2705 | HALYARD* Protectapad, Large | 28.5cm x 43cm | 100 sheets/pack 4 packs/case 400 sheets/case |

Our trusted range of products provide hygienic protection for surfaces and patients with a proven track record of consistent high quality, efficacy and performance.

SURFACES

HALYARD ISOWIPE* Bactericidal Wipes
HALYARD VERSA* Towel
HALYARD TERI* Wipers
HALYARD TERI* Hydro Absorbent Wipers
HALYARD* Absorbent Tray Liner
HALYARD* Benchroll

PATIENTS

HALYARD* Underpads
HALYARD* Bed Protector
HALYARD* Protectapad
HALYARD* Bedsheet
HALYARD* Washcloth Extra
HALYARD* Neutral Hand Wipes

JOYA
MEDICAL SUPPLIES



Advantage*
HALYARD

Customers around the world choose Halyard Health because of our clinically-superior products, purposeful innovation, remarkable service and deep expertise.

It's what we call the HALYARD ADVANTAGE*



HALYARD

O&M Halyard Australia Pty. Limited

Phone: 1800 664 227

www.halyardhealth.com.au

1. Mechanical break resistance test data on file.
Note: These tests are directional based on small sample sizes.

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.

*Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2023 O&M Halyard, Inc. All Rights Reserved.

O&M Halyard Australia Pty Ltd, Suite 1, Level 1, Trinit 3, 39 Delhi Road, North Ryde, NSW 2113, Australia

