



HIGHLY ABSORBENT & VERSATILE PROTECTION

HALYARD* Underpads

TRUSTED PROTECTION

FOR SURFACES & PATIENTS

SO YOU CAN CARE WITH CONFIDENCE

Looking for a product that's hygienic for patients, as well as soft and absorbent for patient wellbeing? A product that will protect all types of surfaces, linen and furniture in a variety of settings such as hospital, aged care, GPs, podiatry, chiropractic, veterinarian and obstetrics? A product that's versatile, waterproof and available in various absorbency levels?

OUR SOLUTION IS HALYARD* UNDERPADS

HALYARD* Underpads are constructed from a soft top layer and highly absorbent layers of tissue with a strong blue polyethylene backing that is quiet, drapable, and helps prevent strikethrough of liquid. Underpads are the ideal choice as a medical procedure underlay where fluids are involved, providing hygienic protection for patients and surfaces. HALYARD* Underpads also assist in reducing the time taken for staff to clean and change linen.

APPLICATIONS

- Protection for linen and mattress
- Bedsheet
- Light incontinence aid
- Chair and surface protection
- Pillow protector
- Underlay for dressings, procedures and surgical preparation
- Underlay for bedpan or cover to prevent splashing when removing
- Baby scale liner

DEPARTMENTS

- Hospitals
- Aged Care Facilities
- Wards
- Obstetrics
- Outpatients
- Casualty
- Children's Wards
- Nursery
- Blood Banks
- GPs
- Allied Health
- Veterinary
- Early Childhood Clinics



HALYARD* Underpads

- ✓ **HIGHLY ABSORBENT¹ TISSUE LAYERS**
Available in 5 ply and 8 ply to suit various absorbency needs
- ✓ **STRONG POLYETHYLENE BACKING**
Helps prevent strikethrough of liquid to keep surfaces dry
- ✓ **SOFT TOP LAYER**
Comfortable, gentle on patients' skin
- ✓ **SOFT, SMOOTH POLY BASE**
Quiet to assist with patient discretion and drapable to protect varied surfaces
- ✓ **BLUE FOLDED EDGES**
Helps to contain fluid by preventing run off, blue edge is a visual guide for edge of pad against white bedsheet
- ✓ **DISPOSABLE**
Hygienic, helps prevent cross-contamination and also assists in reducing the time taken for staff to clean and change linen

JOYA
MEDICAL SUPPLIES



Advantage*
HALYARD

Customers around the world choose Halyard Health because of our clinically-superior products, purposeful innovation, remarkable service and deep expertise.

It's what we call the HALYARD ADVANTAGE*

DISPENSING

Available in two levels; 5 ply and 8 ply, HALYARD* Underpads have a variety of applications including bed sheet / mattress protection, chair protection, pillow protection and bench cover.

- Inner poly bags – make it easier to store on shelves or on trolley, easier to divide up a shipper when smaller quantities are required. Each poly bag will have a label with a barcode, product code, product description and lot number allowing product to be identified easily outside the shipper. See below for packaging details.
- Easier to transport – no need to carry a large shipper, each poly bag can be transported throughout your facility.

PRODUCT INFORMATION

Code	Description	Size	Packaging
2765	HALYARD* Underpads 5 ply	39.5cm x 55cm	4 x 75 sheets per polybag 300 sheets/case
2768	HALYARD* Underpads 8 ply	39.5cm x 55cm	4 x 50 sheets per polybag 200 sheets/case



Our trusted range of products provide hygienic protection for surfaces and patients with a proven track record of consistent high quality, efficacy and performance.

SURFACES

HALYARD ISOWIPE* Bactericidal Wipes
HALYARD VERSA* Towel
HALYARD TERI* Wipers
HALYARD TERI* Hydro Absorbent Wipers
HALYARD* Absorbent Tray Liner
HALYARD* Benchroll

PATIENTS

HALYARD* Underpads
HALYARD* Bed Protector
HALYARD* Protectapad
HALYARD* Bedsheet
HALYARD* Washcloth Extra
HALYARD* Neutral Hand Wipes



O&M Halyard Australia Pty. Limited
Phone: 1800 664 227
www.halyardhealth.com.au

1. Absorbency test data on file.
Note: These tests are directional based on small sample sizes.

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.

*Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2023 O&M Halyard, Inc. All Rights Reserved.

O&M Halyard Australia Pty Ltd, Suite 1, Level 1, Trinita 3, 39 Delhi Road, North Ryde, NSW 2113, Australia