













FitRight Range Guide



Product Selection Guide

Type of Product	Absorbency	Ideal Use	Mobility Level
 Drypads	Maximum	<ul style="list-style-type: none"> For individuals with severe urinary output or wound/fluid management For individuals with high risk for skin breakdown 	
 Bladder Control Pads	Light-moderate	<ul style="list-style-type: none"> For individuals with light to moderate urinary output Can walk with or without assistance 	
 Liners & knit pants	Moderate-maximum	<ul style="list-style-type: none"> For individuals with moderate to severe urinary output and/or bowel incontinence Can walk with or without assistance or may be used for bedridden individuals 	
 Protective Underwear	Light-maximum	<ul style="list-style-type: none"> For individuals with light to severe urinary output with occasional bowel incontinence episodes Can walk with or without assistance 	
 Briefs	Heavy-maximum	<ul style="list-style-type: none"> For individuals with heavy to severe urinary output and/or bowel incontinence Bedridden or limited mobility 	

Mobility Key



Suitable for persons with limited mobility



Suitable for ambulatory persons



Suitable for non ambulatory persons

Types and description	Item No.	Size	Pkg.	Total capacity (ml)
DRYPADS				
Ultrasorbs	ULTRASORBS1016	Ultrasorbs AP 25x41cm	10/bg, 100/CTN	1110
	EXTRASORB2336A	Extrasorb AP 51x91cm	10/bg, 70/CTN	1494
	EXTRASORB3036A	Extrasorb AP 76x91cm	5/bg, 70/CTN	1812
Ultrasorbs Advanced+	ULTRASORBS2436	Ultrasorbs AP 61-91cm	10/bg, 70/CTN	1896
	ULTRASORBS3136	Ultrasorbs AP 79-91cm	10/bg, 40/CTN	2419
Ultrasorbs Advanced Premium with SlipResist Technology	MSC383500	Ultrasorbs AP 76 x 91cm	5/bg, 70/CTN	1840
Medline Skin Fold Sheets	MSC10614	Medline Dry Sheet 15 x 35 cm	10/bg, 250/CTN	256
BLADDER CONTROL PADS				
ContourPlus Bladder Control Pads	BCPE01	Moderate, 14 x 26 cm	28/bg, 336/CTN	281
	BCPE02	Maximum, 16 x 34 cm	28/bg, 168/CTN	465
	BCPE03	Ultimate, 20 x 43 cm	28/bg, 168/CTN	762
FitRight ActivEdge™ Male Guards	MSCMGO2	14 x 25 cm	52/bg, 208/CTN	440
FitRight Double Up Booster Pads	MSC326015	9 x 29 cm	24/bg, 192/CTN	659
FitRight Double Up Booster Pads	MSC90180	18 x 42cm	30/bg, 180/CTN	1894



Medline International Two Australia Pty Ltd.
2 Fairview Place, Marsden Park NSW 2765.
1800 110 511
medline.com.au

Types and description	Item No.	Size	Pkg.	Total capacity (ml)
LINERS & KNIT PANTS				
FitRight Extra Liners	FITLINER100	33 x 76 cm	20/bg, 80/CTN	1554
FitRight Ultra Liners	FITLINER300	33 x 76 cm	20/bg, 80/CTN	1973
FitRight Super Liners	FITLINER500	33 x 76 cm	20/bg, 80/CTN	2363
FitRight Bariatric Liner	ULTRASOFTBARI	52 x 84 cm	18/bg,72/CTN	2392
Knit Pants	MSC86300	Medium/Large, 51-152 cm waist	5/bg, 100/CTN	-
	MSC86400	X-Large, 114-178 cm waist	5/bg, 100/CTN	-
	MSC86500	XX-Large,127-191 cm waist	5/bg, 100/CTN	-
	MSC86700	3XL, 165-216 cm waist	2/bg, 50/CTN	-
PROTECTIVE UNDERWEAR				
FitRight Ultra Underwear	FIT23005A	Medium, 71-102 cm	20/bg, 80/CTN	846
	FIT23505A	Large, 102-142 cm	20/bg, 80/CTN	851
	FIT23600A	X-Large, 142-173 cm	20/bg, 80/CTN	1073
FitRight Super Underwear	FIT33005A	Medium, 71-102 cm	20/bg, 80/CTN	902
	FIT33505A	Large, 102-142 cm	20/bg, 80/CTN	980
	FIT33600A	X-Large, 142-173 cm	20/bg, 80/CTN	1056
FitRight Protection Plus Underwear	MSC33255	Small, 51-71 cm	22/bg, 88/CTN	1060
	MSC33005	Medium, 71-102 cm	20/bg, 80/CTN	1090
	MSC33505	Large, 102-142 cm	18/bg, 72/CTN	1148
	MSC33600	X-Large, 142-173 cm	14/bg, 56/CTN	1301
	MSC33700	XX-Large, 173-203 cm	12/bg, 48/CTN	1107
FitRight Extended Wear™ Underwear	MSC53005	Medium, 71-102 cm	16/bg, 64/CTN	1940
	MSC53505	Large, 102-142 cm	14/bg, 56/CTN	2006
	MSC53600	X-Large, 142-173 cm	12/bg, 48/CTN	2118
FitRight Specialty Size Underwear	FIT700A	XX-Large, 173-203 cm	20/bg, 80/CTN	834
	FIT800	3XL, 191-239 cm	20/bg, 80/CTN	1689
YOUTH & CHILDREN				
Medline Training Pants	HSP29817	Medium, 2-3T, 9-15 kg	17/bg, 136/CTN	668
	HSP29815	Large, 3-4T, 15-19 kg	15/bg, 120/CTN	834
	MSC29813	X-Large, 4-5T, 18+ kg	13/bg, 104/CTN	853
Medline Youth Underwear	MSC23001A	Small/Medium, 8-31 kg	15/bg, 60/CTN	1077
	MSC23003A	Large/X-Large, 27-56 kg	12/bg, 48/CTN	1254
BRIEFS				
FitRight Ultra Briefs	FITULTRAMD	Medium, 81-112 cm	20/bg, 80/CTN	1801
	FITULTRALG	Large, 112-142 cm	20/bg, 80/CTN	1710
	FITULTRAXLG	X-Large, 142-163 cm	20/bg, 80/CTN	1916
	FITULTRAXXL	XXLarge, 152-178 cm	20/bg, 80/CTN	1845
FitRight Stretch Ultra Briefs	FRSU1	Size 1, 76-132 cm	20/bg, 80/CTN	2006
	FRSU2	Size 2, 130-178 cm	20/bg, 80/CTN	2192
FitRight Stretch Extended Wear™ Briefs	FREW1	Size 1, 79-132 cm	20/bg, 80/CTN	2660
	FREW2	Size 2, 102-178 cm	20/bg, 80/CTN	2660
FitRight Extended Wear™ Briefs	MTB80200	Small, 48-79 cm	30/bg, 120/CTN	2552
	MTB80300	Medium, 69-109 cm	15/bg, 60/CTN	3340
	MTB80500	Large, 104-150 cm	15/bg, 60/CTN	4176
	MTB80600	X-Large, 145-165 cm	15/bg, 60/CTN	4345
FitRight Super Specialty Size Briefs	BARIBRIEF	4XL, 165-239 cm	8/bg, 32/CTN	2371
	BARIBRIEFPLUS	5XL, 229-269 cm	8/bg, 32/CTN	2945



Medline International Two Australia Pty Ltd.
 2 Fairview Place, Marsden Park NSW 2765.
 1800 110 511
 medline.com.au