



## Comfort Suction Unit

Model: V3

### Triple Protection consist of

① Overflow valve ② Antibacterial filter ③ Fool-proof design



Patent Nos:  
ZL 2016 3 0298084.1  
Other patents pending

#### It's Comfort

- High performance pump delivers -550mmHg vacuum and 16 liter/min flow rate

#### It's Protected "Triple Protection"

- Overflow valve plus antibacterial filter and fool-proof design that prevents back flow of fluid and airborne contaminants from entering the pump

#### It's Convenient

- Easy to store using tubing holder and accessory compartment

#### It's Ergonomic

- Comprehensive dashboard with ergonomic view angle ensures a smooth operation

#### Technical Specification

- Suction class: High Vacuum / Low Flow
- Max vacuum @ pump inlet: -73 kPa/ -550 mmHg
- Max air flow @ pump inlet: 16 L/min
- Adaptor Input: AC 100-240V, 50/60Hz
- Adaptor Output: DC 12V 2.0A
- Collection canister capacity: 1400ml (Polycarbonate, PC)
- Power Consumption: ≤ 24W
  - Subject to technical modification without prior notice

Model	Qty per carton	Carton volume
V3	4 pcs	0.131 cbm/ctn

#### Tubing holder

##### Tubing Holder



##### Accessory compartment



##### Ergonomic design



Distributed by:

**Silverline Medical Pty Ltd**

Unit 8A, 380 Pennant Hills Road, Pennant Hills, NSW 2120  
1800 959 777

[www.rossmax.com.au](http://www.rossmax.com.au)