



Aneroid Sphygmomanometer

Model : **GB series**
AGC



GB100



- White enamel 300 mmHg no-pin stop manometer
- Chrome-plate brass air release valve
- Latex free Cuff
- Single head Stethoscope



**White enamel 300mmHg
non-stop pin reliable
manometer**



**Chrome-plate brass air
release valve**



**Light-weight
stethoscope**



Latex free cuff

GB102 - D-ring cuff with single head stethoscope



Rossmax Cuff with D-ring :
Adult: 24 - 32cm (9.4" - 12.6")
Other size cuffs available upon request

Model	Qty per carton	Carton volume
GB series	20 pcs	0.042 cbm/ctn

Storage and Transportation Condition
-10°C~60°C(14°F~140°F)
10%~90% RH, 700~1060 hPa

Distributed by:
Silverline Medical Pty Ltd
8A, 380 Pennant Hills Road, Pennant Hills, NSW 2120
1800 959 777
www.rossmax.com.au