

# CONSUMABLES

---



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

THEATRE CONSUMABLES

All surgical consumables are pre-packed and sterile, ready to be used in the operating room.

LATEX-FREE BANDS

ID indicates Inner Diameter. BW indicates Band Width.



DEF1405

DEF1590A	#16 – (ID) 60mm x 1.5mm (BW) Salmon Colour, Packet of 2, Double Packed	100
DEF1793A	#33 – (ID) 90mm x 3mm (BW) Salmon Colour, Packet of 2	25
DEF3542	#33 – (ID) 90mm x 3mm (BW) Salmon Colour, Packet of 2	100
DEF3560	#33 – (ID) 90mm x 3mm (BW) Salmon Colour, Packet of 3	100
DEF669A	#33 – (ID) 90mm x 3mm (BW) Salmon Colour, Packet of 5 Double Packed	25
DEF5797	#33 – (ID) 90mm x 3mm (BW) Blue, Packet of 10	50
DEF619	#34 – (ID) 100mm x 3mm (BW) Clear, Packet of 5	25
DEF600	#34 – (ID) 100mm x 3mm (BW) Clear, Packet of 10	100
DEF624A	#64 – (ID) 90mm x 6mm (BW) Salmon Colour, Packet of 1	100
DEF5528	#64 – (ID) 90mm x 6mm (BW) Clear, Packet of 5	100
DEF626A	#64 – (ID) 90mm x 6mm (BW) Salmon Colour, Packet of 5 Double Packed	100

CLAMPS/CLIPS/HOLDERS

TUBE CLAMPS

DEF1575A	Tube Clamp with Lock Blue Plastic 12cm	50
----------	--	----



DEF1575A

DEF2190	Towel Clamp Green Plastic 11.5cm	25
---------	----------------------------------	----



DEF2190

DEF2488	Towel Clip Stainless Steel 11cm Sharp/Sharp	10
---------	---	----



DEF2488

DEF2298	Tube Clamps 13cm Smooth Jaw	10
---------	-----------------------------	----



DEF2298

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## CLAMPS/CLIPS/HOLDERS (CONTINUED)

### DRAPE CLIP

General purpose, reusable Drape/Peach Clip, used during surgery for clamping purposes to ensure items remain where placed. A wide easy to use grip surface ensures heavy duty holding power.



DEROYAL-DRAPE-CLIP

DEROYAL-DRAPE-CLIP	Assorted Colours Non Sterile	8
--------------------	------------------------------	---




UTILITY-CLIP-UC-1-P, UTILITY-CLIP-UC-1-R

UTILITY-CLIP-UC-1-P	Utility Clip Pink Non Sterile	1
UTILITY-CLIP-UC-1-R	Utility Clip Red Non Sterile	1
DEF6030	Utility Clip Pink Sterile	50
DEF6031	Utility Clip Red Sterile	50

### TUBE HOLDER

Ideal for securing and holding tubes during a procedure. Adhesive backing for securement onto a drape or towel. Adjustable velcro to allow gentle and easy placement.



DEF1889

DEF1889	Tube Holder with Velcro Adhesive 7.5cm x 7.5cm	25
---------	---	----

## COMBINE DRESSING

Highly absorbent layer of cellulose blend enclosed in a soft and conformable non woven fabric. Each non woven combine dressing is pre-cut and each end is sealed into a variety of sizes.



DEF711

DEF684	9cm x 10cm	300
DEF711	9cm x 20cm	25
DEF719	9cm x 20cm, Packet of 5	100
DEF720	9cm x 20cm, Packet of 10	50
DEF697	9cm x 90cm	25
DEF3727	10cm x 15cm	25
DEF3728	10 cm x 20cm	25
DEF3729	20cm x 20cm	15
DEF715	20cm x 20cm	25
DEF80	20cm x 60cm	100
DEF1160	20cm x 80cm	50

## COMBINE ROLL

DEF714	Non Woven Combine Roll 9cm x 10m, Double Packed	16
DEF1004	Non Woven Combine Roll 20cm x 2.5m	25

## COMBINE SPECIALISED



DEF804

DEF804	Non Woven Amputation Dressing 50cm x 70cm	50
--------	--	----

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
<div>DENTAL</div>					
<div>BITE BLOCK</div>					
100% cotton gauze with a barium sulfate detectable thread.					
<div></div>					
DEF1828					
DEF1828	Bite Block 7.5cm X-ray with Tail	100			
DEF1246	Bite Block 7.5cm X-ray with Tail Non Sterile	100			
DEF1232	Bite Block 7.5cm X-ray without Tail Non Sterile	100			
<div>TOOTHBRUSH</div>					
<div></div>					
DEF690					
DEF690	Toothbrush	25			

<div>DENTAL PRESSURE PACK</div>					
<div></div>					
DEF799					
DEF799	X-ray with 30cm Cotton Tail	100			
DEF1063	X-ray with 35cm Tape, Packet of 2	100			
DEF144	3cm x 4cm X-ray Small with 10cm tail, Packet of 2	50			
<div></div>					
DEF144					
DEF145	4cm x 5cm X-ray Large with 9cm tail, Packet of 2	50			
<div></div>					
DEF145					
DEF1209	Dental Roll 4mm x 35mm, Packet of 5	100			
<div></div>					
DEF1209					



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## DEFOGGER

Specially formulated solution that prevents lens fogging. Comes with a non-abrasive, x-ray detectable foam pad. Sterile.



DEROYAL-28-0101

DEROYAL-28-0101	DeFogger Oval Bottle	20
-----------------	----------------------	----

## EYE

### EYE SPEAR

Highly absorbent.



DEF484

DEF484	Eye Spear, Packet of 2	100
DEF485	Eye Spear, Packet of 5 Double Packed	20
DEF486	Eye Spear, Packet of 6 Double Packed	20
DEF487	Eye Spear, Packet of 10 Double Packed	20

### EYE STROLL



DEF112


DEF112	Eye Stroll, Packet of 5	50
DEF081	Eye Stroll, Packet of 10 Double Packed	50

## ELECTROSURGICAL

### DIATHERMY TIP CLEANER

This coated scratch pad cleaner quickly and easily removes debris from instrument tips and diathermys. The adhesive backing will secure the Tip Cleaner onto draped surfaces and contains a radiopaque strip to ensure X-ray detectability.

DEF1886	Tip Cleaner 5cm x 5cm with Adhesive and Radiopaque Strip	50
---------	--	----



DEF1886

### QUIVER

A single use disposable polypropylene Quiver used to hold diathermy during surgical procedures. Quiver is recyclable (metal ring needs to be removed before being placed in recycling container).

DEF2772	5cm x 25cm with Metal Ring	90
---------	----------------------------	----



DEF2772

DEF3068	5cm x 25cm with Cable Tie	90
---------	---------------------------	----



DEF3068

DEF3291	9cm x 35cm with Metal Ring	40
DEF3296	9cm x 35cm with Cable Tie	40
DEF3905	4.5cm x 15.5cm Neuro	100
DEF3906	4.5cm x 20.5cm Neuro	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

### FOAM PAD DRESSING

Polyurethane foam pads for use in a clinical setting.



DEF792

DEF792	100mm x 100mm x 13mm	100
DEF769	100mm x 100mm x 25mm	50
DEF783	100mm x 180mm x 13mm	50
DEF784	100mm x 180mm x 25mm	40
DEF738	150mm x 150mm x 13mm	50



DEF738

DEF1738	150mm x 315mm x 25mm Non Sterile	50
DEF914	150mm x 200mm x 25mm	25
DEF1549	200mm x 25mm x 25mm	100



DEF1549

DEF913	200mm x 300mm x 25mm	25
DEF739	250mm x 250mm x 13mm	50
DEF786	300mm x 600mm x 13mm Double Packed	20

### PREP APPLICATORS

DEF921	Prep Foam 50mm x 60mm x 25mm White No Handle	100
--------	---	-----



DEF921

DEF1846	Prep Foam on a Handle 15cm Blue	100
---------	------------------------------------	-----



DEF1846

DEF1847	Prep Foam on a Handle 15cm Blue, Packet of 2	100
DEF2035	Prep Foam on a Handle 20cm Orange	100



DEF2035

### PREP SETS

DEF3106	PREP SET 2 x Gallipot 150ml 3 x Kidney Dish 700ml Clear 1 x Kidney Dish 700ml Yellow Wrapped in a 91cm x 91cm Field	20
DEF3458A	PREP SET 5 x Gauze Swabs 10cm x 10cm XRD 2 x Sponge Holder 19cm Plastic 2 x Gallipot 150ml 1 x Kidney Dish 700ml Clear Wrapped in a 91cm x 115cm Field	30

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## IDENTILOOPS

These radio-opaque Identiloops are used in surgical procedures for retraction, identification or occlusion of arteries, veins, nerves or tendons. A soft and flexible, inert medical grade silicone which is oval in cross-section, allows for even pressure distribution to minimise damage to tissues. Available in 4 colours – Blue, White, Red and Yellow in mini and maxi sizes. All double packed.

2200220-2	Mini-Red 40cm, Packet of 2	10
2200220-4	Mini-Red 40cm, Packet of 4	10



2200220

2200420-2	Mini-Yellow 40cm, Packet of 2	10
-----------	-------------------------------	----



2200420

2200620-2	Mini-Blue 40cm, Packet of 2	10
2200620-4	Mini-Blue 40cm, Packet of 4	10



2200620

2200920-2	Mini-White 40cm, Packet of 2	10
2200920-4	Mini-White 40cm, Packet of 4	10



2200920

2202220-2	Maxi-Red 40cm, Packet of 2	10
2202220-4	Maxi-Red 40cm, Packet of 4	10



2202220

2202420-2	Maxi-Yellow 40cm, Packet of 2	10
2202420-4	Maxi-Yellow 40cm, Packet of 4	10



2202420

2202620-2	Maxi-Blue 40cm, Packet of 2	10
-----------	-----------------------------	----



2202620

2202920-2	Maxi-White 40cm, Packet of 2	10
-----------	------------------------------	----



2202920

DEF1193	Super Maxi-Blue 40cm, Packet of 2	10
DEF1318	Mixed Maxi 40cm, Packet of 5	10
DEF2182	Super Maxi-Blue 60cm, Packet of 2	10




DEF2182

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------


### LIGHT HANDLE COVERS

Our range of flexible light handle covers are single sterile and come either coloured as blue or green. These adaptable covers are twisted or clicked onto the light handle adapter and are firmly held for duration required. We can supply the compatible light handle adapters for each operating room light. Easy and versatile to use.




DEF2082

DEF2082	Light Handle Cover Green Friction Fit	50
DEF2919	Light Handle Cover Green Friction Fit, Packet of 2	40
DEF2675	Light Handle Cover Blue, Click On	50




DEF2675

DEF3170	Light Handle Cover Blue Click On, Packet of 2	40
DEF2557	Light Handle Cover Foam Blue	50



DEF2557


DEF3673	Camera Light Handle Cover 75mm	25
---------	--------------------------------	----



DEF3673


### IRRIGATION BULB SYRINGE

AS011P	Irrigation Bulb Syringe 60ml	50
--------	------------------------------	----




AS011P

### SKIN HOOKS



DEF336A

DEF336A	5 Small Hooks in a Needle Case Double Packed	20
DEF1032A	5 Large Hooks in a Needle Case Double Packed	20




DEF2402

DEF2402	Catspaw Retractor, Packet of 2	10
---------	--------------------------------	----

### YARD STRIPS

DEF858	1/2 Yard Strips (53cm) X-ray Packet of 5	50
--------	--	----



DEF858



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## NEEDLE COUNTERS DEROYAL


Excellent selection of Needle Counters from DeRoyal with easy to read numbers and description, facilitating correct counting of blades and needles used in a procedure. Full surface magnetic face provides extra space for sharps, pins, blade or needles.

'Y' codes indicate Needle Counter case is Yellow.

### DELUXE NEEDLE COUNTERS

#### MAGNET

25-0103Y	Single Magnet, Clear Top 10 Count, 20 Capacity	50
----------	---	----



25-0103Y

#### MAGNET/FOAM COMBINATION

25-0401	Magnet/Foam 20 Count, 40 Capacity	20
---------	--------------------------------------	----




25-0401

25-0701	Foam Strips/Magnet 15 Count, 30 Capacity	50
---------	---	----



25-0701

25-0702	Magnet/Foam Strips 30 Count, 60 Capacity	27
---------	---	----



25-0702

#### FOAM

25-0504	Foam/Adhesive 40 Count, 80 Capacity	12
---------	--	----



25-0504

25-0803	Foam Strips 60 Count, 120 Capacity	27
---------	---------------------------------------	----











25-0803

## TOUCH FASTENER

DEF5826	5cm x 16cm Velcro Strip X-ray Detectable, Pack of 2	100
---------	--	-----



DEF5826

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.																																																												
<div>DEBLADE NEEDLE COUNTERS</div> <div>This range of Needle Counters all feature a blade removal system which provides simple removal of contaminated blades.</div> <div>MAGNET</div> <table><tr><td>25-1104Y</td><td>Double Magnet with DeBlade 30 Count, 60 Capacity</td><td>50</td><td colspan="3"></td></tr></table> <div></div> <div>25-1104Y</div> <table><tr><td>25-1108</td><td>Double Magnet with DeBlade 60 Count, 120 Capacity</td><td>27</td><td colspan="3"></td></tr></table> <div></div> <div>25-1108</div> <div>MAGNET/FOAM COMBINATION</div> <table><tr><td>25-1401</td><td>Magnet/Foam Block with DeBlade 20 Count, 40 Capacity</td><td>20</td><td colspan="3"></td></tr><tr><td>25-1401Y</td><td>Magnet/Foam Block with DeBlade 20 Count, 40 Capacity</td><td>20</td><td colspan="3"></td></tr></table> <div></div> <div>25-1401</div> <table><tr><td>25-1701</td><td>Magnet/Foam Strips with DeBlade 15 Count, 30 Capacity</td><td>50</td><td colspan="3"></td></tr></table> <div></div> <div>25-1701</div> <table><tr><td>25-1702</td><td>Magnet/Foam Strips with DeBlade 30 Count, 60 Capacity</td><td>27</td><td colspan="3"></td></tr></table> <div></div> <div>25-1702</div> <div>MARKING PENCIL</div> <table><tr><td>DEF464</td><td>Prepacked Marking Pencil – Omnichrom 108-9</td><td>24</td><td colspan="3"></td></tr></table> <div></div> <div>DEF464</div> <div>SAFETY PINS</div> <table><tr><td>DEF071SP</td><td>#2 (38mm)</td><td>50</td><td colspan="3"></td></tr><tr><td>DEF635</td><td>#3 (51mm) Bulk Non Sterile</td><td>144</td><td colspan="3"></td></tr></table> <div></div> <div>DEF635</div> <table><tr><td>DEF073SP</td><td>#4 (57mm)</td><td>50</td><td colspan="3"></td></tr></table> <div></div> <div>DEF073SP</div>						25-1104Y	Double Magnet with DeBlade 30 Count, 60 Capacity	50				25-1108	Double Magnet with DeBlade 60 Count, 120 Capacity	27				25-1401	Magnet/Foam Block with DeBlade 20 Count, 40 Capacity	20				25-1401Y	Magnet/Foam Block with DeBlade 20 Count, 40 Capacity	20				25-1701	Magnet/Foam Strips with DeBlade 15 Count, 30 Capacity	50				25-1702	Magnet/Foam Strips with DeBlade 30 Count, 60 Capacity	27				DEF464	Prepacked Marking Pencil – Omnichrom 108-9	24				DEF071SP	#2 (38mm)	50				DEF635	#3 (51mm) Bulk Non Sterile	144				DEF073SP	#4 (57mm)	50			
25-1104Y	Double Magnet with DeBlade 30 Count, 60 Capacity	50																																																															
25-1108	Double Magnet with DeBlade 60 Count, 120 Capacity	27																																																															
25-1401	Magnet/Foam Block with DeBlade 20 Count, 40 Capacity	20																																																															
25-1401Y	Magnet/Foam Block with DeBlade 20 Count, 40 Capacity	20																																																															
25-1701	Magnet/Foam Strips with DeBlade 15 Count, 30 Capacity	50																																																															
25-1702	Magnet/Foam Strips with DeBlade 30 Count, 60 Capacity	27																																																															
DEF464	Prepacked Marking Pencil – Omnichrom 108-9	24																																																															
DEF071SP	#2 (38mm)	50																																																															
DEF635	#3 (51mm) Bulk Non Sterile	144																																																															
DEF073SP	#4 (57mm)	50																																																															

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## SKIN MARKING PENS

Skin Markers are used during surgical procedures and contain medical grade, non irritating gentian violet ink. They mark a smooth, clear line on a patient's skin, and the ink is non toxic, non tattooing and quick drying.

26-001	Skin Marker Standard Tip with 15cm Ruler	50
--------	--	----



26-001

26-002	Skin Marker Fine Tip with 15cm Ruler	50
--------	--------------------------------------	----



26-002

26-001PR	Skin Marker Standard Tip Prep Resistant Ink with 15cm Ruler	50
----------	---	----

26-003PR	Skin Marker Standard Tip Prep Resistant Ink with 15cm Ruler and 8 labels (55mm x 11mm each)	50
----------	---	----

26-005PR	Skin Marker Standard Tip Prep Resistant Ink	50
----------	---	----



26-001PR

## NON-STERILE

26-018NS	Mini Pre-Op Skin Marker Bulk Non Sterile	200
----------	--	-----



26-018NS

## UTILITY MARKER

Ideal for writing on instruments, specimens, vials and tubes, basins, bowls, and other objects as needed within the sterile field. This marker is not to be used on skin.

DEF3219	Utility Marker Fine Tip Black Ink	25
---------	-----------------------------------	----



DEF3219

## SPONGES & SWABS

### SPONGES ABDOMINAL X-RAY

Soft highly absorbent cotton sponges. Each complete with an x-ray radiopaque element for easy detection.



DEF162

DEF162	23cm x 23cm with Tape, XRD	100
DEF160	23cm x 23cm with Tape, XRD Packet of 3	100
DEF168	23cm x 23cm with Tape, XRD Packet of 5	100
DEF1351	30cm x 30cm, XRD, Packet of 5 Double Packed	75
DEF1352	30cm x 48cm, XRD, Packet of 5 Double Packed	60



DEF1351

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## GAUZE SWABS X-RAY

Excellent absorbency and quality swabs with 8 ply cotton blend. Includes an x-ray radiopaque element for easy detection.



DEF044

DEF202	7.5cm x 7.5cm X-ray, Packet of 5 Double Packed	100
DEF044	10cm x 10cm X-ray, Packet of 5 Double Packed	100
DEF997	10cm x 10cm X-ray, Packet of 5 32ply, Double Packed	50

## GAUZE SWABS GREEN



DEF190

DEF190	10cm x 10cm Green, Packet of 2	100
DEF192	10cm x 10cm Green, Packet of 3	100
DEF207	10cm x 10cm Green, Packet of 5 Double Packed	100
DEF225	10cm x 10cm Green, Packet of 20	50
DEF196	10cm x 10cm Green, Bulk Non Sterile	100

## GAUZE SWABS FOLDED

DEF347	Gauze Swab 7cm x 21cm X-ray, Folded Packet of 5	50
--------	---	----



DEF347

DEF185	Gauze Swab Neuro 10cm x 10cm 32ply X-ray, Packet of 10, Double Packed	100
--------	---	-----



DEF185

DEF1074	Gauze Swab Neuro 10cm x 20cm X-ray Folded, Packet of 5	25
---------	--	----

## GAUZE HEAD DRESSINGS

DEF611	Head Dressing Kit includes: 8 x Gauze Swabs 10cm x 20cm 20 x Gauze Swabs 10cm x 10cm 2 x 90cm x 90cm Gauze, Folded	35
--------	---	----

## GAUZE PREP SPONGES X-RAY

Round sponges with X-ray radiopaque element. Excellent absorbency. Used as a prep sponge prior to a surgical procedure.



DEF251

DEF255	Sponges Prep X-ray, 6cm x 6cm	100
DEF251	Sponges Prep X-ray, 6cm x 6cm Packet of 2	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## SPONGES & SWABS (CONTINUED)

### GAUZE PEANUT SPONGES

Round sponges ideal for delicate surgical sponging or soft tissue dissection. X-ray radiopaque element for easy detection.

30-106	5 in a Foam Holder, XRD	50
--------	-------------------------	----



30-106

### SPECIALISED GAUZE



DEF169

DEF986	Donor Site Gauze Dressing, 5cm x 10cm, Double Packed	100
DEF169	Recovery Gauze 5cm x 5cm Plain ,32 Ply	100

## SCOPE WARMER

Long, double sided Scope Warmer from DeRoyal that pre-warms endoscopes to body temperature. Easy to use design provides 1 hour of warming time on each side.

DEROYAL-TM500D	Scope Warmer	10
----------------	--------------	----



DEROYAL-TM500D

## PVC TUBING WITH XRD

Medical grade, clear and transparent tubing. A radiopaque strip along length of tubing ensures X-ray detectability.

DEF3051	PVC Tubing with XRD 16Fr x 17cm, Packet of 6 Double Packed	25
---------	--	----



DEF3051

DEF3052	PVC Tubing with XRD 16Fr x 22cm, Packet of 6 Double Packed	25
---------	--	----

## SCALPEL BLADE REMOVER

Single handed blade remover that can hold up to three used surgical blades. Simple to use with a quick 'push, click and pull' action to remove and contain used scalpel blades.

QSSVCAS-3Y	Qlicksmart Blade Cassette	100
------------	---------------------------	-----



QSSVCAS-3Y

Please note - blade not included.



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

SURGICAL RETRACTOR

ONETRAC

ONETRAC is a single-use cordless surgical retractor with a built-in LED light source and smoke evacuation channel. It provides illumination in the depths of a surgical tissue pocket or cavity. No additional parts or components or fibre optic cables are required. Ideal for breast or reconstructive procedures, as well as other surgical procedures requiring access, illumination and retraction of soft tissue. Reduce cost, increase efficiency and eliminate unnecessary risk of cross-contamination with ONETRAC.



ONETRAC

C100150-5	40mm x 20mm	5
C100170-5	55mm x 8mm	5
C100140-5	90mm x 22mm	5
C100180-5	105mm x 60mm	5
C100165-5	135mm x 20mm	5
C100110-5	135mm x 30mm	5
C100190-5	175mm x 30mm	5



ONETRAC LX

Features integrated multi-LED light sources and dual smoke evacuation channels.

C100270	30mm x 15mm & 50mm x 20mm Dual Blade	1
C100240	90mm x 20mm	1
C100270	90mm x 30mm	1
C100230	90mm x 60mm	1
C100290	135mm x 20mm	1
C100210	135mm x 30mm	1
C100260	135mm x 60mm	1
C100220	175mm x 30mm	1



C100260, C100230

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## ESMARCH BANDAGE

Offers excellent compression and elastic properties to achieve exsanguination, whilst delivering protection for latex sensitive patients and health professionals.

DEF1006	6.5cm x 2m Green Silicone Double Packed	10
---------	--	----



DEF1006

DEF2764	10cm x 2.7m Blue Neoprene Double Packed	20
---------	--	----



DEF2764


35-309	7.6cm x 2.7m Blue Neoprene	20
35-409	10cm x 2.7m Blue Neoprene	20
35-609	15cm x 2.7m Blue Neoprene	20



35-309

## SUTURE BOOTS WITH ADHESIVE TRAY

Suture boots offer a secure way to hold sutures without slippage or damage during a surgical procedure. They are fitted onto the distal working end of surgical instruments to provide an improved atraumatic grip and even pressure distribution when clamping an artery or vein. The tray comes with an adhesive strip for securing and easy removal. Radio-opaque element ensures X-ray detection. Available in a Yellow colour.



30-702

DEF1093	2 on an Adhesive Tray Double Packed	25
DEF1094	4 on an Adhesive Tray Double Packed	25
30-702	6 on an Adhesive Tray	50
DEF1096	8 on an Adhesive Tray Double Packed	25
30-701	10 on an Adhesive Tray	50
DEF2284	12 on an Adhesive Tray Double Packed	25

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## SILICONE TUBING

Our 100% medical grade, clear silicone tubing range is available pre-cut and sterile or as non-sterile rolls. They are double packed unless otherwise specified.

Sizing below is represented as:

ID (Internal Diameter) x OD (Outer Diameter) x Length of tube in cm or m.

DEF2307	0.30mm x 0.64mm x 20cm	75
DEF2522	1mm x 2mm x 50cm	50
DEF477	1.57mm X 2.41mm X 30cm Packet of 2	25
DEF1545	1.57mm x 3.18mm x 1m	15
DEF1428	1.57 mm x 2.41mm x 15cm	100
DEF473	1mm x 3mm x 50cm, Packet of 4	20
DEF3649	1.98mm x 3.18mm x 75cm Packet of 2	20
DEF5712-10	2mm x 4mm x 45cm	10
DEF2314	3.18mm x 6.35mm x 25cm	60
DEF2038A-10	4mm x 6mm x 45cm	10
DEF1665	5mm x 7mm x 20cm Packet of 2	25
DEF2167B-10	5mm x 7mm x 45cm	10
DEF1253	5mm x 8mm x 1.6m	50
DEF1261	5mm x 9mm x 1.5m Single Packed	50
DEF5655-10	6mm x 9mm x 45cm	10
DEF2786	6.4mm x 9.5mm x 2cm Packet of 2	375
DEF2386	6.4mm x 9.5mm x 1m	15
DEF5652	6.4mm x 11mm x 3m	8
DEF1818A-10	7mm x 10mm x 45cm	10
DEF2032	7mm x 10mm x 3m	8
DEF482	7mm x 12mm x 85cm Packet of 2	100
DEF1859A-10	8mm x 12mm x 45cm	10
DEF1204	8mm x 14mm x 3m	5
DEF2468	10mm x 16mm x 1.5m	16

## SILICONE TUBING WITH CONNECTOR

Flexible pre-cut and sterile silicone tubing with connector attached. All of this range is double packed.



DEF1186

DEF735	5mm x 8mm x 1m	50
DEF1186	5mm x 9mm x 30cm	50
DEF306A	5mm x 9mm x 50cm	50
DEF1712	5mm x 9mm x 60cm	40
DEF1337	5mm x 9mm x 1m	50
DEF1231A	5mm x 9mm x 1.5m	30
DEF305A	7mm x 11mm x 1m	30

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## WOUND DRAINAGE

Wound drains and evacuators provide drainage for various surgical procedures, allowing for the removal of blood and excess fluid to improve patient outcome. Our wound drainage systems are made up of evacuators, flat drains, fluted drains, round drains, capillary action drains and easy flow drains.

They are all made from 100% medical grade silicone, which is inert, soft and flexible and designed to reduce the potential of post – operative infections.

### WOUND EVACUATOR BULB

Our wound evacuator is a clear bulb with easy to read graduations up to 100ml/cc, and is designed to slowly draw fluid and blood from the wound site. It creates an ideal vacuum to a maximum of 125mmHg pressure and is compatible with our flat drains, fluted drains and round drains.



31300380DS

31300380DS	Evacuator Bulb 100ml/cc	10
------------	-------------------------	----

### PERFORATED FLAT WOUND DRAIN

Soft and comfortable silicone flat drains with multiple drainage holes to enhance the fluid removal in a wound. They have an ergonomic flat shape which enables an easier placement. X-ray detectable.



11100200780DS

11100200780DS	7mm x 800mm with a 220mm Channel, White XRD Strip	10
11100100780DS	7mm x 800mm with a 220mm Channel, Blue XRD Strip	10
11100201080DS	10mm x 800mm with a 220mm Channel, White XRD Strip	10
11100101080DS	10mm x 800mm with a 220mm Channel, Blue XRD Strip	10


### FLUTED FLAT WOUND DRAIN

Clear soft silicone drain with white x-ray detectable end.

11101400780DS	7mm x 1020m with 200mm Channel and Trocar	10
11101411080DS	10mm x 1020mm with 200mm Channel and Trocar	10

## FLUTED ROUND WOUND DRAIN

Fluted drains are made of soft, flexible silicone with or without attached trocar. Designed with a multi channel to allow for great drainage and assisting in minimal tissue in-growth. The inner core of the H design offers drain strength and flexibility. X-ray detectable.



11003021080DS	10Fr x 1100mm with 300mm Channel and Connector, Blue XRD Strip	10
11003121080DS	10Fr x 1100mm with 200mm Channel and Connector, Blue XRD Strip	10
11003121580DS	15Fr x 1200mm with 300mm Channel and Connector, Blue XRD Strip	10
11003121980DS	19Fr x 1200mm with 300mm Channel and Connector, Blue XRD Strip	10
11003122480DS	24Fr x 1100mm with 300mm Channel and Connector, Blue XRD Strip	10

### FLUTED ROUND WOUND DRAIN SETS

11003131080DS	10Fr x 1100mm with 200mm Channel and Trocar, Blue XRD Strip	10
11003131580DS	15Fr x 1200mm with 300mm Channel and Trocar, Blue XRD Strip	10
11003131980DS	19Fr x 1200mm with 300mm Channel and Trocar, Blue XRD Strip	10
11003831080DS	10Fr x 1100mm with 300mm Channel and Trocar, White XRD Strip	10
11003811580DS	15Fr x 1100mm with 300mm Channel and Trocar, White XRD Strip	10
11003811980DS	19Fr x 1100mm with 300mm Channel and Trocar, White XRD Strip	10

### PERFORATED ROUND WOUND DRAIN

Clear soft silicone drain with a blue x-ray detectable line.

11000301280DS	12Fr 500mm with 5 holes 2x4mm	10
11000301880DS	18Fr 500mm with 5 holes 4x6mm	10
11000302480DS	24Fr 500mm with 5 holes 5x8mm	10
11000302780DS	27Fr 500mm with 5 holes 6x9mm	10
11000303080DS	30Fr 500mm with 5 holes 7x10mm	10
11000303680DS	36Fr 500mm with 5 holes 8x12mm	10
11000600780DS	7Fr 1250mm with 72 Holes	10
11000601080DS	10Fr 1250mm with 72 Holes	10
11000601580DS	15Fr 1250mm with 72 Holes	10
11000601980DS	19Fr 1250mm with 72 Holes	10

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## MEDICAL CONSUMABLES

### ANAESTHETIC FACE MASK

Our anaesthetic face masks are ergonomically designed to optimise the application of anaesthetics. The masks have an inflating air cushion that is adjustable by an inflation valve, for sealing and patient comfort.

The transparent face plate and cushion provides easy patient facial visibility. The masks are for single patient use and are provided non-sterile.



A1701

A1701	Size 0 - Neonate	50
A3011	Size 0 - Neonate, Strawberry Scented	50
A1702	Size 1 - Infant	50
A3012	Size 1 - Infant, Strawberry Scented	50
A1703	Size 2 - Toddler	50
A3013	Size 2 - Toddler, Strawberry Scented	50
A1704	Size 3 - Child / Small Adult	50
A3014	Size 3 - Child / Small Adult Strawberry Scented	50
A1705	Size 4 - Medium Adult	50
A3015	Size 4 - Medium Adult, Strawberry Scented	50
A1706	Size 5 - Large Adult	50
A3016	Size 5 - Large Adult, Strawberry Scented	50
A1707	Size 6 - X-Large Adult	50
A3017	Size 6 - X-Large Adult, Strawberry Scented	50

### BURNS

#### BURNS DRESSING

Specialised products for patients with burns. Non adherent, highly absorbent wound dressings.



DEF987

DEF199	10cm x 20cm, Packet of 20	50
DEF679	50cm x 70cm, Melolin Double Packed	35
DEF3570	50cm x 175cm Melolin Double Packed	25
DEF970	90cm x 50cm	100
DEF974	90cm x 2m	50



DEF679

#### BURNS GAUZE

Soft roll suitable for wounds and burns.

DEF965	10cm x 80cm, Packet of 16	50
DEF2436	25cm x 90cm, Packet of 2	350
DEF1039	60cm x 90cm, Packet of 6	100
DEF342	90cm x 2m	50
DEF344	90cm x 2m, Packet of 4	25
DEF1662	90cm x 5m	10

#### BURNS PLASTIC BAG

DEF650	305mm x 390mm	100
--------	---------------	-----



CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## CLIP SEAL BAG



DEF3286

DEF3286	20cm x 25cm, Packet of 20	50
DEF1703	28cm x 38cm	25

## FREEZER BAG



DEF639

DEF639	30cm x 45cm	100
DEF66	1 x 30cm x 45cm Bag with 2 x Latex-free Bands #34 Double Packed	100
DEF1336	2 x 30cm x 45cm Bags with 2 x Latex-free Bands #34 Double Packed	100

## DRAWSTRING BAG



DEF846

DEF846	28cm x 54.5cm, Double Packed	25
DEF864	44.5cm x 46.5cm	25
DEF861	44.5cm x 46.5cm	100
DEF1687	44.5cm x 46.5cm, Packet of 2	50
DEF848	50cm x 50cm	25
DEF847	50cm x 85cm	25

## PLASTIC BAG



DEF1878

DEF1878	14cm x 30cm	100
DEF1276	14cm x 30cm with a Blank Adhesive Label	100
DEF647	20cm x 20cm	100
DEF638	28cm x 32.5cm	100
DEF650	Burns 30.5cm x 39cm	100
DEF642	45.5cm x 51cm with Cotton Tape 60cm Long Double Packed	100
DEF808	Heavy Duty 57.5cm x 62cm	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## COTTON PRODUCTS

### COTTON BALLS

Packaged and sterilised for ease of entry into a sterile field. Ideal for prepping, wound cleaning, perineal care and fluid application.



DEF074

DEF074	3cm, Packet of 5	100
--------	------------------	-----

### COTTON WOOL ROLL

Soft and absorbent cotton wool roll used for padding in patient care. Sterilised and ready to be used.



DEF874

DEF874	15cm x 55cm, Double Packed	50
DEF875	30cm x 90cm, Double Packed	25
DEF880	38cm x 165cm, Double Packed	15

## SPONGES

### SPONGES NASAL BOLSTER

Used after nasal or sinus surgery to absorb any discharge or ooze from the nose. High quality non woven net filled with absorbent cotton padding, and stretch elastic sides to secure around the ears. Extremely soft and comfortable for post surgery patients. Latex free and available as sterile or non sterile.



DEF1381

DEF1381	2cm x 8cm with Ear Loops	20
DEF1380	2cm x 8cm with Ear Loops Bulk Non Sterile	20
DEF1382	2cm x 8cm with Ear Loops	100

## NON WOVEN PRODUCTS

### NON WOVEN BALLS

DEF1354	2.5cm, Packet of 5	100
---------	--------------------	-----

### NON WOVEN GAUZE SWABS - STERILE

DEF3933A	5cm x 5cm, Packet of 2	100
DEF1131	5cm x 5cm, Packet of 5	100

### NON WOVEN GAUZE SWABS - NON STERILE

DEF1869	5cm x 5cm, Bulk	100
DEF1868	7.5cm x 7.5cm, Bulk	100
DEF1867	10cm x 10cm, Bulk	100

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## SWABS

### GAUZE SWABS PLAIN - STERILE

Sterile and ready to be used for a wide variety of applications, including as a dressing, for prepping procedures and wound cleaning. Available as a sterile Plain 8 ply gauze.



DEF181

DEF183	7.5cm x 7.5cm, Packet of 3	100
DEF181	7.5cm x 7.5cm, Packet of 5	100
DEF371	7.5cm x 7.5cm, Packet of 20	50
DEF189	10cm x 10cm, Packet of 2	100
DEF186	10cm x 10cm, Packet of 3	100



DEF045

DEF045	10cm x 10cm, Packet of 5	100
DEF208	10cm x 10cm, Packet of 10	50
DEF213	10cm x 10cm, Packet of 10	300
DEF178	10cm x 10cm, Packet of 20	25
DEF177	10cm x 10cm, Packet of 25	25
DEF214	10cm x 10cm, Packet of 25	200
DEF1033	10cm x 20cm, Packet of 4	100
DEF198	10cm x 20cm, Packet of 5	250
DEF195	10cm x 20cm, Packet of 10	50
DEF853	10cm x 80cm, Packet of 4 16 ply Folded	50

### GAUZE SWABS PLAIN – NON STERILE

DEF180	5cm x 5cm, Bulk	100
DEF197	7.5 x 7.5cm, Bulk	100
DEF136	10cm x 10cm, Bulk	100

### GAUZE DRESSING

DEF571	15cm x 22cm	100
--------	-------------	-----

### PURCOTTON SWABS

Made from a non-woven fabric using natural cotton, our purcotton swabs have a high liquid absorbency and are 100% biodegradable. They are also low linting and have an excellent wet-strength.

DEF5753	10cm x 10cm, Pack of 10	100
DEF1062	10cm x 10cm, Non Sterile	100

## SANITARY DISPOSAL BAG

DEF709	Packet of 100, Non Sterile	20
--------	----------------------------	----



DEF709

## VAGINAL PLUG



DEF850SP

DEF850SP	90mm x 30mm x 30mm	100
----------	--------------------	-----


CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

### SWAB STICKS/COTTON TIPPED APPLICATORS


Packaged and sterilised, ready to use for taking a swab or specimen from a patient. Paper, plastic or wooden shaft in different length sizes with a single or double ended soft cotton head.

#### SWAB STICKS WITH WOODEN SHAFT



DEF057

DEF057	7.5cm Single Ended, Packet of 3	100
DEF442	7.5cm Single Ended, Packet of 5	100
DEF439	15cm Single Ended, Packet of 2	100
DEF059	15cm Single Ended, Packet of 3	100
DEF443	15cm Single Ended, Packet of 5	100




DEF439

#### SWAB STICKS WITH WOODEN SHAFT – NON STERILE

DEF062	15cm Single Ended, Bulk	1000
--------	-------------------------	------

#### SWAB STICKS WITH PAPER SHAFT




DEF421

DEF058	7.5cm Double Ended, Packet of 3	100
DEF421	7.5cm Double Ended, Packet of 10	100


### SWAB STICKS WITH PAPER SHAFT – NON STERILE

DEF061	7.5cm Double Ended, Bulk	1250
--------	--------------------------	------




DEF061

#### SWAB STICKS JUMBO



DEF401

DEF400	15cm, Packet of 2	150
DEF401	15cm, Packet of 3	150
DEF402	15cm, Packet of 5	150
DEF1593	20cm	200
DEF455	20cm, Packet of 2	150
DEF407	20cm, Packet of 3	150



DEF407

#### SWAB STICKS JUMBO – NON STERILE

DEF055	15cm Bulk	200
DEF054	20cm Bulk	200

CODE	DESCRIPTION	BOX QTY.	CODE	DESCRIPTION	BOX QTY.
------	-------------	----------	------	-------------	----------

## SILICONE TUBING ROLL – NON STERILE

Excellent range of medical grade silicone tubing available in a variety of different metric sizes. Available as non sterile rolls of 15 or 25 metres.



2110030239	0.30mm (ID) x 0.64mm (OD) 15m Roll
2110050134	0.5mm (ID) x 1mm (OD) 25m Roll
203005009	0.51mm (ID) x 0.94mm (OD) 15m Roll
203006011	0.64mm (ID) x 1.19mm (OD) 15m Roll
2110070439	0.7mm (ID) x 1.6mm (OD) 15m Roll
201010020	1mm (ID) x 2mm (OD) 25m Roll
2110100234	1mm (ID) x 3mm (OD) 25m Roll
2110100439	1.016mm (ID) x 2.159mm (OD) 15m Roll
203014019	1.47mm (ID) x 1.96mm (OD) 15m Roll
2110150434	1.5mm (ID) x 3mm (OD) 25m Roll
2110150939	1.57mm (ID) x 2.41mm (OD) 15m Roll
2110151039	1.57mm (ID) x 3.18mm (OD) 15m Roll
2110190239	1.98mm (ID) x 3.18mm (OD) 15m Roll
201020030	2mm (ID) x 3mm (OD) 25m Roll
2110200334	2mm (ID) x 4mm (OD) 25m Roll
2110200534	2mm (ID) x 6mm (OD) 25m Roll
2110260239	2.65mm (ID) x 4.88mm (OD) 15m Roll
2110300234	3mm (ID) x 5mm (OD) 25m Roll
2110300334	3mm (ID) x 6mm (OD) 25m Roll
2110300434	3mm (ID) x 7mm (OD) 25m Roll
2110310239	3.18mm (ID) x 6.35mm (OD) 15m Roll
2110330239	3.35mm (ID) x 4.65mm (OD) 15m Roll

2110400434	4mm (ID) x 6mm (OD) 25m Roll
2110400634	4mm (ID) x 7mm (OD) 25m Roll
2110400734	4mm (ID) x 8mm (OD) 25m Roll
2110470139	4.76mm (ID) x 7.93mm (OD) 15m Roll
20304707901	4.77mm (ID) x 7.9mm (OD) 15m Roll
20304709501	4.77mm (ID) x 9.50 (OD) 15m Roll
2110500234	5mm (ID) x 7mm (OD) 25m Roll
2110500334	5mm (ID) x 8mm (OD) 25m Roll
2110500534	5mm (ID) x 9mm (OD) 25m Roll
2110500734	5mm (ID) x 10mm (OD) 25m Roll
2110600234	6mm (ID) x 9mm (OD) 25m Roll
2110600334	6mm (ID) x 10mm (OD) 25m Roll
2110600434	6mm (ID) x 11mm (OD) 25m Roll
2110600534	6mm (ID) x 12mm (OD) 25m Roll
2110630139	6.4mm (ID) x 9.5mm (OD) 15m Roll
2110630239	6.4mm (ID) x 11mm (OD) 15m Roll
2110630339	6.4mm (ID) x 12.7mm (OD) 15m Roll
2110700334	7mm (ID) x 10mm (OD) 25m Roll
2110700534	7mm (ID) x 11mm (OD) 25m Roll
2110700634	7mm (ID) x 12mm (OD) 25m Roll
201070130	7mm (ID) x 13mm (OD) 25m Roll
2110790239	7.9mm (ID) x 12.7mm (OD) 15m Roll
2110800334	8mm (ID) x 12mm (OD) 25m Roll
2110800534	8mm (ID) x 14mm (OD) 25m Roll
2110900334	9mm (ID) x 13mm (OD) 25m Roll
2110950139	9.5mm (ID) x 12.7mm (OD) 15m Roll
2110950239	9.5mm (ID) x 14.3mm (OD) 15m Roll
2110950339	9.5mm (ID) x 15.9mm (OD) 15m Roll
2111000334	10mm (ID) x 14mm (OD) 25m Roll
2111000434	10mm (ID) x 15mm (OD) 25m Roll
201100160	10mm (ID) x 16mm (OD) 25m Roll
201100180	10mm (ID) x 18mm (OD) 25m Roll
2111200334	12mm (ID) x 17mm (OD) 25m Roll
2111250134	12.5mm (ID) x 19mm (OD) 25m Roll
2111270239	12.7mm (ID) x 19.05mm (OD) 15m Roll