



Adaptors and Drain Valves

Urocare® Replacement Adaptors, Valves and Gasket-Rings

The following contain all replacement adaptors, valves and gasket-rings which are included in Urocare® products that are sold separately. If you are unable to find an item, chances are the item is listed under another category or it is not sold separately.

Urocare® Male Urinal to Latex Leg Bag Adaptor

This replacement adaptor allows for a direct connection between Urocare® Reusable Latex Leg Bags and Urocare® Male Urinal Sheaths. Only for use with Urocare® products.

Sold by the each, 1 adaptor per package.

This is a stocked item.

Non-Sterile, Latex-Free

LATEX-FREE



LEARN MORE



Item No.	Description	Dimensions†	Case Item No.
6022	Standard Top Adaptor	1.24" x 0.84" (3.15cm x 2.13cm)	602212 (12/CS)

Urocare® Straight-Thru Adaptor with Thumb Clamp

This short, multi-purpose adaptor may be used as a replacement adaptor for the Standard Bottom Drain Valve (#6004) supplied with all Urocare® Reusable Latex Urinary Leg Bags and Urocare® Male Urinal Sheaths.

It may serve as an alternative drain valve, an extension between Urocare® Male Urinal Sheaths and Urocare® reusable latex urinary leg bags or Uro-Safe® disposable vinyl urinary leg bags or as an adaptor for Urocare® Male Urinal Sheaths to connect them directly to most any reusable or disposable urinary leg bag, urinary drainage bag or other urinary and ostomy appliances which accept tubing as a method for connection. Includes small, black Gasket-Ring (#5999) and Two-Position Thumb Clamp (#7000) preassembled.

Sold by the each, 1 adaptor per package.

This is a stocked item.

Non-Sterile, DEHP & Latex-Free.

LATEX-FREE



LEARN MORE



Item No.	Description	Dimensions†	Case Item No.
6005	Straight-Thru Adaptor with Clamp	0.31" I.D. x 5.5" (8mm Ø x 14cm)	600512 (12/CS)

Urocare® Standard Bottom Drain Valve

Features sealed threads so fingers don't come in contact with urine while draining. Only requires 1½ turns to open or close. Replacement drain valve for use with all Urocare® Reusable Latex Leg Bags and Urocare® Male Urinal Sheaths. Includes Small Black Gasket-Ring (#5999) preassembled.

Sold by the each, 1 valve per package.

This is a stocked item.

Non-Sterile, Latex-Free.

LATEX-FREE



LEARN MORE



Item No.	Description	Dimensions†	Case Item No.
6004	Standard Bottom Drain Valve	1.53" x 0.6" (3.88cm x 1.52cm)	600410 (100/CS)

Urocare® Urinary Drainage Bottle Adaptor

Provides a secure connection for Urocare® Clear-Vinyl Drainage Tubing. Replacement adaptor for use with Urocare® Urinary Drainage Bottle. May also be used with Urocare® Reusable Latex Leg Bags as a replacement for the Standard Top Adaptor (#6002) for those who use clear-vinyl tubing.

Sold by the each, 1 adaptor per package.

This is a stocked item.

Non-Sterile, Latex-Free.

LATEX-FREE



LEARN MORE



Item No.	Description	Dimensions†	Case Item No.
6003	Drainage Bottle Adaptor	1.24" x 0.84" (3.15cm x 2.13cm)	600310 (10/CS)

† Approximate dimensions in inches(in) and centimeters(cm).