

SUPRAPUBIC CATHETER



FROM SAYCO PTY LTD.

- ◆ Shortened at the tip to reduce potential trauma
- ◆ Open ended tip +2 eyelets to maximise drainage
- ◆ Catheter Placement Markings
- ◆ In-Line balloon reduces ridging and cuffing



Unit 5, 260 Wickham Rd, Moorabbin, Victoria 3189
Ph: 03 9553 2700 or 1800 333 456 Fax: 03 9553 4858
Email: customerservice@sayco.net.au Web: www.sayco.com.au
SupraPlus® is a registered trademark of Sayco® Pty Ltd.

Supra Plus®

- ◆ Sizes 14Fr~24Fr, Silicone Catheter
- ◆ Catheter Placement Markings for easier changeover
- ◆ 41cm length
- ◆ 10ml balloon capacity (14Fr is 5ml)
- ◆ Radio opaque line
- ◆ Supra-pubic and Urethral use



sayco

Specifications

Size	Length	Quantity per box	Code	Class IIb	Open Ended Tip	Shortened Tip	Balloon Capacity	In-Line Balloon
14FR	41cm	10	SP14FR	✓	✓	✓	5ml	✓
16FR	41cm	10	SP16FR	✓	✓	✓	10 ml	✓
18FR	41cm	10	SP18FR	✓	✓	✓	10 ml	✓
20FR	41cm	10	SP20FR	✓	✓	✓	10 ml	✓
22FR	41cm	10	SP22FR	✓	✓	✓	10 ml	✓
24FR	41cm	10	SP24FR	✓	✓	✓	10 ml	✓

Supra Plus® is a high quality catheter that can be used as a SupraPubic and Urethral Catheter